

24162	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24163	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24164	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24165	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24166	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24167	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24168	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24169	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit .0
24170	Replacement Frits	for Passive Air Sampling Kit 3pk
24171	Replacement Frits	for Passive Air Sampling Kit 3pk, Siltek Treated
24172	TO-Can Canister	1 Liter, 1/4" Valve
24173	TO-Can Canister	3 Liter, 1/4" Valve (cat #563647 is needed if connecting to a sampling kit)
24174	TO-Can Canister	6 Liter, 1/4" Valve
24174-650	TO-Can Canister	6 Liter, 1/4" Valve w/Silcosteel
24174-653	TO-Can Canister	6 Liter, 1/4" Valve Ship Under Vacuum
24175	TO-Can Canister	15 Liter, 1/4" Valve
24176	TO-Can Canister	1 Liter w/Gauge, 1/4" Valve
24176-650	TO-Can Canister	1 Liter w/Gauge, 1/4" Silcosteel Valve
24177	TO-Can Canister	3 Liter w/Gauge, 1/4" Valve (cat #563647 is needed if connecting to a sampling kit)
24178	TO-Can Canister	6 Liter w/Gauge, 1/4" Valve
24178-652	TO-Can Canister	6 Liter 1/4" Valve 2" McDaniels Vacuum Pressure Gauge 30" HG to 30 PSI
24179	TO-Can Canister	15 Liter w/Gauge, 1/4" Valve
24180	SilcoCan Canister	1 Liter, 1/4" Valve
24180-650	SilcoCan Canister	1 Liter w/Siltek Treated 1/4" Valve
24181	SilcoCan Canister	3 Liter, 1/4" Valve (cat #563646 is needed if connecting to a sampling kit)
24181-650	SilcoCan Canister	3 Liter w/Siltek Treated 1/4" Valve (cat #563646 is needed if connecting to a samp
24182	SilcoCan Canister	6 Liter, 1/4" Valve
24182-650	SilcoCan Canister	6 Liter w/Siltek Treated 1/4" Valve
24182-65065	SilcoCan Canister	6 Liter w/Siltek 1/4" Valve ship under vacuum
24183	SilcoCan Canister	15 Liter, 1/4" Valve
24183-650	SilcoCan Canister	15 Liter w/Siltek Treated 1/4" Valve
24183-65065	SilcoCan Canister	15 Liter, 1/4" Silcosteel Valve ship under vacuum
24185	Fittings	Quick Connect Stem Fitting for the miniature Air Sampling Canister
24186	Fittings	Quick Connect Stem Fitting Sulfinert Treated, for the Miniature Air Sampling Canis
24187	Fittings	Quick Connect Body Fitting for the Miniature Air Sampling Canister
24188	Miniature Canister	400cc Electro Polished With Quick Connect Stem Fitting
24189	Miniature Canister	400cc Sulfinert With Quick Connect Stem Fitting
24190	Miniature Canister	400cc Sulfinert With Sulfinert Quick Connect Stem Fitting
24191	Miniature Canister	400cc Electro Polished With Metal Seated Diaphragm Valve
24192	Miniature Canister	400cc Sulfinert With Metal Seated Diaphragm Valve
24193	Miniature Canister	400cc Sulfinert With Sulfinert Treated Diaphragm Valve
24194	Miniature Canister	1000cc Electro Polished With Quick Connect Stem Fitting
24195	Miniature Canister	1000cc Sulfinert With Quick Connect Stem Fitting
24196	Miniature Canister	1000cc Sulfinert With Sulfinert Quick Connect Stem Fitting
24197	Miniature Canister	1000cc Electro Polished With Metal Seated Diaphragm Valve
24198	Miniature Canister	1000cc Sulfinert With Metal Seated Diaphragm Valve
24199	Miniature Canister	1000cc Sulfinert With Sulfinert Treated Diaphragm Valve
24205	Miniature Canister	400cc Electro Polished w/Nut and Ferrule
24206	Miniature Canister	1000cc Electro Polished w/Nut and Ferrule
24207	Miniature Canister	400cc Sulfinert w/Nut and Ferrule
24208	Miniature Canister	1000cc Sulfinert w/Nut and Ferrule
24215	Canister Carrying Box Kit	4 x 6-Liter Carrying Boxes w/Plastic Handles, 1 x 4pk Shipper
24216	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Non-Silcosteel, Complete Kit, .(
24217	Sampling Kit	CONFIRM DELIVERY WITH OPERATIONS > 25 Siltek, Complete Kit .0008 Orific
24218	Orifice	Replacement, .5-2sccm Non-Siltek
24219	Orifice	Replacement, .5-2sccm Siltek Treated
24226	Deluxe Carrying Case	for Two 6 Liter Canisters (14lbs, 20x14x14") Ships One per Box
24229	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflow 423XL, 0.5-2sccm orific
24230	Ultra Clean Resin	100 gm bottle
24230	Ultra Clean Resin	100 gm bottle
24230	Ultra Clean Resin	100 gm bottle
24232	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflow 423XL, 0.5-2sccm orifi
24233	Orifice	Replacement, 2-4 sccm Siltek-treated
24234	Orifice	Replacement, 4-8 sccm Siltek-treated
24235	Orifice	Replacement, 8-15 sccm Siltek-treated
24236	Orifice	Replacement, 15-30 sccm Siltek-treated

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034