

24237	Orifice	Replacement, 30-80 sccm Siltek-treat
24238	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 2-4sccm orif. Silt
24239	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 4-8sccm orif. Silt
24245	Orifice	Replacement, 2-4 sccm Non-Siltek
24246	Orifice	Replacement, 4-8 sccm Non-Siltek
24247	Orifice	Replacement, 8-15 sccm Non-Siltek
24248	Orifice	Replacement, 15-30 sccm Non-Siltek
24249	Orifice	Replacement, 30-80 sccm Non-Siltek
24255	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 2-4sccm orif. Silt
24256	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 4-8sccm orif. Silt
24257	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 8-15 sccm orif. S
24258	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 15-30 sccm orif.
24259	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 30-80 sccm orif.
24260	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 2-4sccm orif. No
24261	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 4-8sccm orif. No
24262	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 8-15 sccm orif. N
24263	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflow 423XL, 15-30 sccm ori
24264	Flow Controller	CONFIRM DELIVERY WITH OPERATIONS > 25 Veriflo 423XL, 30-80 sccm orif.
24265	Siltek 7um Filter	7um In-Line Filter, Siltek Treated(5um no longer manufactured)
24266	Filter	7um In-Line Filter, Non-Silcosteel (5um no longer manufactured)
24267	Air Sampling Timer	for Air Canisters
24268	Test Gauge	Digital, 3" Diameter, 1/4"NPT, 30"Hg to 60psig
24269	Gauge	30" to 0 psig Vacuum 1/8" NPT CBM, 2" diam, SS
24270	Gauge	30" to 0 psig Vacuum 1/4" NPT CBM, 2" diam, SS
24282	Humidification Chamber	
24285	Test Gauge	Analog, 6" Diameter, 1/4" NPT, 30"Hg to 60psig
24295	PUFs	Precleaned PUF plug for ORBO-2000 Assembly
24298	Rupture Disc	Replacement Rupture Disc 2850psi Substitute 26325
24323	EPC Test Kit	Used for EPC Leak Checking. For Use w/Agilent 6890 GCs
24324	Deactivated Wool	10 g
24351	Vials, Crimp Seal	11mm Blue PTFE/RR Pack of 100
24352	Vials, Crimp Seal	11mm Blue PTFE/RR Pack of 1000
24353	Vials, Crimp Seal	11mm Green PTFE/RR Pack of 100
24354	Vials, Crimp Seal	11mm Green PTFE/RR Pack of 1000
24355	Vials, Crimp Seal	11mm Red PTFE/RR Pack of 100
24356	Vials, Crimp Seal	11mm Red PTFE/RR Pack of 1000
24357	Vials, Crimp Seal	11mm Yellow PTFE/RR Pack of 100
24358	Vials, Crimp Seal	11mm Yellow PTFE/RR Pack of 1000
24359	Aluminum Crimp Seals	w/Septa, 11mm Silver PTFE/SIL Pack of 100
24360	Aluminum Crimp Seals	w/Septa, 11mm Silver PTFE/SIL Pack of 1000
24361	Vials, Crimp Seal	11mm Blue PTFE/SIL Pack of 100
24362	Vials, Crimp Seal	11mm Blue PTFE/SIL Pack of 1000
24363	Vials, Crimp Seal	11mm Green PTFE/SIL Pack of 100
24364	Vials, Crimp Seal	11mm Green PTFE/SIL Pack of 1000
24365	Vials, Crimp Seal	11mm Red PTFE/SIL Pack of 100
24366	Vials, Crimp Seal	11mm Red PTFE/SIL Pack of 1000
24367	Vials, Crimp Seal	11mm Yellow PTFE/SIL Pack of 100
24368	Vials, Crimp Seal	11mm Yellow PTFE/SIL Pack of 1000
24369	Aluminum Crimp Seals	11mm Silver, PTFE/Silicone/PTFE Pack of 100
24370	Aluminum Crimp Seals	11mm Silver, PTFE/Silicone/PTFE Pack of 1000
24383	Vials, Crimp Top	12x32mm/2.0mL/11mm Clr w/Mark Pack of 100
24383-220	Vials, Snap Seal	12x32mm/2.0mL/11mm Clr w/Mark Deact 100PK
24384	Vials, Crimp Top	12x32mm/2.0mL/11mm Clr w/Mark Pack of 1000
24384B	Vials, Snap Seal	12x32mm/2.0mL/11mm Blue w/Mark Pack of 1000
24384G	Vials, Snap Seal	12x32mm/2.0mL/11mm Green w/Mark Pack of 1000
24384R	Vials, Snap Seal	12x32mm/2.0mL/11mm Rust w/Mark Pack of 1000
24384Y	Vials, Snap Seal	12x32mm/2.0mL/11mm Yellow w/Mark Pack of 1000
24385	Vials, Crimp Top	12x32mm/2.0mL/11mm Ambr w/Mark Pack of 100
24386	Vials, Crimp Top	12x32mm/2.0mL/11mm Ambr w/Mark Pack of 1000
24386-221	Vials, Snap Seal	12x32mm/2.0ml/11mm Amber w/mark, 1000pk, Deactivated
24433	Seals, Poly Crimp	Clear, PTFE/Butyl Rubber Liner Pack of 100
24434	Seals, Poly Crimp	Clear, PTFE/Butyl Rubber Liner Pack of 1000
24443	Seals, Poly Crimp	Clear PTFE/Silicone Liner Pack of 100
24444	Seals, Poly Crimp	Clear PTFE/Silicone Liner Pack of 1000
24445	Seals, Poly Crimp	Blue PTFE/Silicone Liner Pack of 100
24446	Seals, Poly Crimp	Blue PTFE/Silicone Liner Pack of 1000

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034