

24453 Seals, Poly Crimp  
 24454 Seals, Poly Crimp  
 24455 Seals, Poly Crimp  
 24456 Seals, Poly Crimp  
 24473 Caps, Short Screw  
 24474 Caps, Short Screw  
 24485 Caps, Short Screw  
 24486 Caps, Short Screw  
 24487 Caps, Short Screw  
 24488 Caps, Short Screw  
 24493 Caps, Short Screw  
 24494 Caps, Short Screw  
 24495 Caps, Short Screw  
 24496 Caps, Short Screw  
 24497 Caps, Short Screw  
 24498 Caps, Short Screw  
 24499 Caps, Short Screw  
 24500 Caps, Short Screw  
 24501 Caps, Short Screw  
 24502 Caps, Short Screw  
 24503 Caps, Short Screw  
 24504 Caps, Short Screw  
 24505 Caps, Short Screw  
 24506 Caps, Short Screw  
 24508 Vial Inserts  
 24510 Vial Inserts  
 24511 Vial Inserts  
 24512 Vial Inserts  
 24513 Vial Inserts  
 24516 Vial Inserts  
 24516-220 Vial Inserts  
 24517 Vial Inserts  
 24517-221 Vial Inserts  
 24518 Vial Inserts  
 24518A Vial Inserts  
 24519 Vial Inserts  
 24520 Vial Inserts  
 24528 Syringe, Hamilton  
 24529 Syringe, Hamilton  
 24530 Syringe, Hamilton  
 24531 Syringe, Hamilton  
 24532 Syringe, Hamilton  
 24533 Syringe, Hamilton  
 24534 Syringe, Hamilton  
 24535 Syringe, Hamilton  
 24536 Syringe, Hamilton  
 24537 Syringe, Hamilton  
 24538 Syringe, Hamilton  
 24539 Syringe, Hamilton  
 24540 Syringe, Hamilton  
 24541 Syringe, Hamilton  
 24542 Syringe, Hamilton  
 24543 Syringe, Hamilton  
 24544 Syringe, Hamilton  
 24545 Syringe, Hamilton  
 24546 Syringe, Hamilton  
 24547 Syringe, Hamilton  
 24548 Syringe, Hamilton  
 24549 Syringe, Hamilton  
 24550 Syringe, Hamilton  
 24551 Syringe, Hamilton  
 24553 Syringe, Hamilton  
 24554 Syringe, Hamilton  
 24555 Syringe, Hamilton  
 24556 Syringe, Hamilton  
 24557 Syringe, Hamilton

Chromalytic Technology Pty Ltd

Page 103  
 Clear PTFE/Sil/PTFE Liner Pack of 1l  
 Clear PTFE/Sil/PTFE Liner Pack of 1l  
 Blue PTFE/Sil/PTFE Liner Pack of 10L  
 Blue PTFE/Sil/PTFE Liner Pack of 1000  
 9mm Blue PTFE/Red Rubber Liner Pack of 100  
 9mm Blue PTFE/Red Rubber Liner Pack of 1000  
 9mm Blue PTFE/Sil Pack of 100  
 9mm Blue PTFE/Sil Pack of 1000  
 9mm Green PTFE/Sil Pack of 100  
 9mm Green PTFE/Sil Pack of 1000  
 9mm Yellow PTFE/Sil Pack of 100  
 9mm Yellow PTFE/Sil Pack of 1000  
 9mm Black PTFE/Sil/PTFE Pack of 100  
 9mm Black PTFE/Sil/PTFE Pack of 1000  
 9mm Blue PTFE/Sil/PTFE Pack of 100  
 9mm Blue PTFE/Sil/PTFE Pack of 1000  
 9mm Green PTFE/Sil/PTFE Pack of 100  
 9mm Green PTFE/Sil/PTFE Pack of 1000  
 9mm Natural PTFE/Sil/PTFE Pack of 100  
 9mm Natural PTFE/Sil/PTFE Pack of 1000  
 9mm Red PTFE/Sil/PTFE Pack of 100  
 9mm Red PTFE/Sil/PTFE Pack of 1000  
 9mm Yellow PTFE/Sil/PTFE Pack of 100  
 9mm Yellow PTFE/Sil/PTFE Pack of 1000  
 100ul Glass LV w/PP Bot Spring Pack of 100  
 For SO Screw & Crimp Glass 250ul Flat Bottom 1000Pack  
 For SO Screw & Crimp Poly 100ul No Spring Pack of 100  
 For SO Screw & Crimp Poly 100ul w/Bot Spring 100Pack  
 For LO Crimp & Seal Glass 50ul w/Bot Spring 100Pack  
 For LO Crimp Top Vial 250ul Glass Flange 100Pack  
 For LO Crimp Top Vial Deact. 250ul Glass Flange 100pk  
 For LO Crimp & Seal Glass 350ul Flat Bottom 1000Pack  
 for LO Crimp & Seal Glass, 1000 Pack Deactivated 350ul Flat Bottom  
 250ul Polypropylene with Bottom Spring 100PK  
 250ul, Graduated Markings, Polypropylene Bottom Spring, 100pk  
 For LO Crimp & Seal Poly. 250ul Plas. Flange 100Pack  
 For LO Crimp & Seal Poly. 250ul No Spring Pack of 100  
 1702 WISP Syringe, 25ul HPLC Autosampler  
 1702 WISP Syringe, 250ul HPLC Autosampler  
 701RN/10ul/26s/2in/2pt  
 702N/25ul/22s/2in/2pt  
 702RN/25ul/22s/2in/2pt  
 705N/50ul/22s/2in/2pt  
 705RN/50ul/22s/2in/2pt  
 710N/100ul/22s/2in/2pt  
 710RN/100ul/22s/2in/2pt  
 725N/250ul/22s/2in/2pt  
 725RN/250ul/22s/2in/2pt  
 750N/500ul/22/2in/2pt  
 750RN/500ul/22/2in/2pt  
 95N/5ul/26s/2in/2pt  
 95RN/R/5ul/26s/2in/2pt  
 901N/10ul/26s/2in/2pt  
 901RN/R/10ul/26s/2in/2pt  
 7000.5KH/0.5ul/25/2.75in/2pt  
 7000.5KH/0.5ul/25/2.75in/3pt  
 7001KH/1ul/25s/2.75in/2pt  
 7001KH/1ul/25s/2.75in/3pt  
 7101KH/1ul/22s/2.75in/2pt  
 7101KH/1ul/22s/2.75in/3pt  
 7002KH/2ul/25/2.75in/2pt  
 7102KH/2ul/23/2.75in/2pt  
 7102KH/2ul/23/2.75in/3pt  
 7105KH/5ul/24/2.75in/2pt  
 7105KH/5ul/24/2.75in/3pt  
 1701N/10ul/26s/2in/2pt/GT

AUSTRALIAN Distributors

03 9762 2034



**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd

**Australian Distributors**  
 Importers & Manufacturers  
 www.chromtech.net.au

Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

