

24558	Syringe, Hamilton	1701RN/100ul/26s/2in/2pt/GT
24559	Syringe, Hamilton	1702N/25ul/22s/2in/2pt/GT
24560	Syringe, Hamilton	1702RN/25ul/22s/2in/2pt/GT
24561	Syringe, Hamilton	1705N/50ul/22s/2in/2pt/GT
24562	Syringe, Hamilton	1705RN/50ul/22s/2in/2pt/GT
24563	Syringe, Hamilton	1710N/100ul/22s/2in/2pt/GT
24564	Syringe, Hamilton	1710RN/100ul/22s/2in/2pt/GT
24565	Syringe, Hamilton	1710LT/100ul/GT
24566	Syringe, Hamilton	1710TLL/100ul/GT
24567	Syringe, Hamilton	1725N/250ul/22s/2in/2pt/GT
24568	Syringe, Hamilton	1725RN/250ul/22/2in/2pt/GT
24569	Syringe, Hamilton	1725LT/250ul/GT
24570	Syringe, Hamilton	1725TLL/250ul/GT
24571	Syringe, Hamilton	1750LTN/500ul/22/2in/2pt/GT
24572	Syringe, Hamilton	1750RN/500ul/22/2in/2pt/GT
24573	Syringe, Hamilton	1750LT/500ul/GT
24574	Syringe, Hamilton	1750TLL/500ul/GT
24575	Syringe, Hamilton	1001LTN/1ml/22/2in/2pt/GT
24576	Syringe, Hamilton	1001RN/1ml/22/2in/2pt/GT
24577	Syringe, Hamilton	1001LT/1ml/GT
24578	Syringe, Hamilton	1001TLL/1ml/GT
24579	Syringe, Hamilton	1001.25LTN/1.25ml/22/2in/2pt/GT
24580	Syringe, Hamilton	1001.25LT/1.25ml/GT
24581	Syringe, Hamilton	1002LTN/2.5ml/22/2in/2pt/GT
24582	Syringe, Hamilton	1002RN 2.5ml GT (22/2"/2pt)
24583	Syringe, Hamilton	1002LT/2.5ml/GT
24584	Syringe, Hamilton	1002TLL/2.5ml/GT
24585	Syringe, Hamilton	1005LTN/5ml/22/2in/2pt/GT
24586	Syringe, Hamilton	1005RN 5ml GT (22/2"/2)
24587	Syringe, Hamilton	1005LT/5ml/GT
24588	Syringe, Hamilton	1010LTN/10ml/22/2in/2pt/GT
24589	Syringe, Hamilton	1010RN 10ml GT(22/2"/2)
24590	Syringe, Hamilton	1010LT/10ml/GT
24591	Syringe, Hamilton	1025TLL/25ml/GT
24592	Syringe, Hamilton	HP AS 75ASN/5ul/26s
24593	Syringe, Hamilton	HP AS 75ASN/5ul/23s-26s
24594	Syringe, Hamilton	HP AS 75ASN/5ul/23s-26s Pack of 6
24595	Syringe, Hamilton	HP AS 701ASN/10ul/26s
24596	Syringe, Hamilton	HP AS 701ASN/10ul/23s-26s
24597	Syringe, Hamilton	HP AS 701ASRN/10ul/26s
24598	Syringe, Hamilton	HP AS 701ASRN/10ul/23s-26s
24599	Syringe, Hamilton	HP AS 701ASN/10ul/26s Pack of 6
24600	Syringe, Hamilton	HP AS 701ASN/10ul/23s-26s Pack of 6
24605	Syr Needles, Hamilton	LT/TLL Metal Hub/25s/1.97"/3pt Pack of 6
24606	Syr Needles, Hamilton	LT/TLL Metal Hub/22s/2in/2pt Pack of 6
24607	Syr Needles, Hamilton	LT/TLL Kel-F/26s/2in/2pt Pack of 6
24608	Syr Needles, Hamilton	LT/TLL Kel-F/25s/2in/3pt Pack of 6
24609	Syr Needles, Hamilton	LT/TLL Kel-F/22s/2in/2pt Pack of 6
24613	Chaney Adapter	For 700/1700/1000/5-10ul
24615	Syringe Guide	For 7000/.5-5ul or 700/1700/1000/25-500ul
24616	Syringe Guide	For 700/1700/1000/5-10ul
24617	Syringe, Hamilton	75RN/5ul/26s/2in/2pt
24619	Vial, Screw Thread	2ml Clr w/White Grad Mk Spot Pack of 100
24620	Vial, Screw Thread	2ml Clr w/White Grad Mk Spot Pack of 1000
24620-221	Vial, Screw Thread	2ml Clr w/White Grad Mk Spot IP Deact., Pack of 1000
24621	Vial, Screw Thread	2ml Ambr w/White Grad Mk Spot Pack of 100
24622	Vial, Screw Thread	2ml Ambr w/White Grad Mk Spot Pack of 1000
24623	Vial, Limited Volume	100ul Poly. w/8mm Thread Pack of 100
24624	Vial, Limited Volume	100ul Poly. w/8mm Thread Pack of 1000
24625	Vial, Limited Volume	100ul Glass/Clr Plastic w/8mm Thread Pack of 100
24626	Vial, Step	2ml Clr w/White Grad Mk Spot Pack of 100
24626-220	Vial, Step	2ml Clr w/White Grad Mk Spot Deactivated, Pk of 100
24626-222	Vial, Step	2ml Clr w/White Grad Mk Spot Siltek Deact. Pack of 100
24627	Vial, Step	2ml Clr w/White Grad Mk Spot Pack of 1000
24628	Vial, Step	2ml Ambr w/White Grad Mk Spot Pack of 100
24628-220	Vial, Step	2ml Ambr w/white Grad Mk Spot Pack of 100, Deactivated

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

 <b>HROMalytic</b> +61(0)3 9762 2034 <b>ECH</b> nology Pty Ltd	<b>Australian Distributors</b> Importers & Manufacturers www.chromtech.net.au
Website NEW : <a href="http://www.chromalytic.net.au">www.chromalytic.net.au</a> E-mail : <a href="mailto:info@chromtech.net.au">info@chromtech.net.au</a> Tel: 03 9762 2034 . . . in AUSTRALIA	

[Ctrl F]  
> Cat.#s

> Current System AUD Prices

