

24703	Syringe, SGE	Model 10R 10ul/26/2
24704	Syringe, SGE	Model 25F 25ul/25/2
24705	Syringe, SGE	Model 25R 25ul/25/2
24706	Syringe, SGE	Model 50F 50ul/25/2
24707	Syringe, SGE	Model 50R 50ul/25/2
24708	Syringe, SGE	Model 100F 100ul/25/2
24709	Syringe, SGE	Model 100R 100ul/25/2
24710	Syringe, SGE	Model 250F 250ul/25/2
24711	Syringe, SGE	Model 250R 250ul/25/2
24712	Syringe, SGE	Model 500F 500ul/25/2
24713	Syringe, SGE	Model 500R 500ul/25/2
24715	Syringe, SGE	Model SK-10F 10ul/26/2 Pack of 6
24716	Syringe, SGE	10ul RN Syr w/Plunger Protect. (6pk) Pack of 6
24717	Syringe, Hamilton	1050TLL 50.0ml Gastight
24718	Syringe, SGE	SuperfleX 10ul/26/FN Pack of 6
24719	Syringe, SGE	SuperfleX 10ul/26/RN Pack of 6
24720	Syringe, SGE	SuperfleX Plunger Model 5FX 5ul/26
24721	Syringe, SGE	SuperfleX Plunger Model 5RX 5ul/26
24722	Syringe, SGE	SuperfleX Plunger Model 10FX 10ul/26
24723	Syringe, SGE	SuperfleX Plunger Model 10RX 10ul/26
24724	Syringe, Hamilton	1100TLL 100.0ml Gastight
24725	Syringe, SGE	Model 10F-GT 10ul/26/2
24726	Syringe, SGE	Model 10R-GT 10ul/26/2
24727	Syringe, SGE	Model 25F-GT 25ul/25/2
24728	Syringe, SGE	Model 25R-GT 25ul/25/2
24729	Syringe, SGE	Model 50F-GT 50ul/25/2
24730	Syringe, SGE	Model 50R-GT 50ul/25/2
24731	Syringe, SGE	Model 50F-LT-GT 50ul
24732	Syringe, SGE	Model 50F-LL-GT 50ul
24734	Syringe, SGE	Model 100F-GT 100ul/25/2
24735	Syringe, SGE	Model 100R-GT 100ul/25/2
24736	Syringe, SGE	Model 00F-LT-GT 100ul
24737	Syringe, SGE	Model 100F-LL-GT 100ul
24738	Syringe, SGE	Model 250F-GT 250ul/25/2
24739	Syringe, SGE	Model 250R-GT 250ul/25/2
24740	Syringe, SGE	Model 50F-LT-GT 250ul
24741	Syringe, SGE	Model 250F-LL-GT 250ul
24742	Syringe, SGE	Model 500F-GT 500ul/25/2
24743	Syringe, SGE	Model 500R-GT 500ul/25/2
24744	Syringe, SGE	Model 00F-LT-GT 500ul
24745	Syringe, SGE	Model 500F-LL-GT 500ul
24746	Syr Needles, SGE	Model NP0.5B-S-0.47 .5ul/26
24747	Syr Needles, SGE	Model NP0.5B-S-0.63 .5ul/23
24748	Syr Needles, SGE	Model N10-S-0.47 10ul/26 Pack of 2
24749	Syr Needles, SGE	Model N10-S-0.63 10ul/23 Pack of 2
24750	Syringe, SGE	Model 1MR-GT 1ml/23/2
24751	Syringe, SGE	Model 1MF-LT-GT 1ml
24752	Syringe, SGE	Model 1MDF-LL-GT 1ml
24753	Syringe, SGE	Model 2.5MDR-GT 2.5ml/23/2
24754	Syringe, SGE	Model 2.5MF-LT-GT 2.5ml
24755	Syringe, SGE	Model 2.5MDF-LL-GT 2.5ml
24756	Syringe, SGE	Model 5MDR-GT 5ml/23/2
24757	Syringe, SGE	Model 5MDR-LL-GT 5ml
24758	Syringe, SGE	Model 10MDR-GT 10ml/23/2
24759	Syringe, SGE	Model 10MDR-LL-GT 10ml
24760	Syringe, SGE	Model 25MDR-LL-GT 25ml
24761	Syringe, SGE	Model 50MR-LL-GT 50ml
24762	Syringe, SGE	Model 100MR-LL-GT 100ml
24763	Syr Needles, SGE	Luer Lock for GasTight Syr Pack of 5
24764	Syr Needles, SGE	Luer Lock for GasTight Syr Pack of 5
24765	Syr Needles, SGE	Luer Lock for GasTight Syr Pack of 2
24766	Syringe, SGE	Guided Plunger Micro Syr Model 5R-GP 5ul/26/2
24767	Syringe, SGE	Guided Plunger Micro Syr Model 10F-GP 10ul/26/2
24768	Syringe, SGE	Guided Plunger Micro Syr Model 10R-GP 10ul/26/2
24769	Syringe, SGE	Guided Plunger Syringe 5ul Model 5F-GP
24770	Syringe, SGE	Model 0.5BR-OC-7/0.47 .5ul/26/1

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

 HROMalytic +61(0)3 9762 2034 ECHnology Pty Ltd		Australian Distributors Importers & Manufacturers www.chromtech.net.au
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA		

[Ctrl F]
> Cat.#s

> Current System AUD Prices

