

24771	Syringe, SGE	Model 0.5BR-7BV .5ul/23/2
24772	Syringe, SGE	Model 0.5BR-7 .5ul/23/1
24773	Syringe, SGE	Model 1BR-5 1.0ul/23/1
24774	Syringe, SGE	Model 1BR-7 1.0ul/23/1
24775	Syringe, SGE	Model 1BR-7BV 1.0ul/23/2
24776	Syringe, SGE	Model 1BR-OC-7/0.47 1.0ul/26/1
24777	Syringe, SGE	Model 5BR-5 5.0ul/23/1
24778	Syringe, SGE	Model 5BR-7 5.0ul/23/1
24779	Syringe, SGE	Model 5BR-7BV 5.0ul/23/2
24780	Syringe, SGE	HP Autosampler Model 5F-HP-0.47 5ul/26/fixed
24781	Syringe, SGE	HP Autosampler Model 5F-HP-0.63 5ul/23/fixed
24782	Syringe, SGE	HP Autosampler, Pack of 6 SK-5F-HP-0.47 5ul/26/fixed
24783	Syringe, SGE	HP Autosampler Pack of 6 ModelSK-5F-HP-0.63 5ul/23fixed
24784	Syringe, SGE	HP Autosampler Model 10F-HP-0.47/10ul/26fixed
24785	Syringe, SGE	HP Autosampler Model 10F-HP-0.63 10ul/23fixed
24786	Syringe, SGE	HP Autosampler, Pack of 6 Model SK-10F-HP-0.47 10ul/26/F
24787	Syringe, SGE	HP Autosampler Pack of 6 Model SK-10F-HP-0.63 10ul/23/F
24789	Syringe, SGE	HP Autosampler Model 10F-HP-GT-0.63 10ul/23/F
24790	Syringe, SGE	HP Autosampler Model .5BR-HP-0.47 .5ul/26/R
24791	Syringe, SGE	HP Autosampler Model .5BR-HP-0.63 .5ul/23/R
24792	Syringe, SGE	HP Autosampler Model 5R-HP-0.47 5ul/26/R
24793	Syringe, SGE	HP Autosampler Model 5R-HP-0.63 5ul/23/R
24794	Syringe, SGE	HP Autosampler Model 10R-HP-0.47 10ul/26/R
24795	Syringe, SGE	HP Autosampler Model 10R-HP-0.63 10ul/23/R
24796	Syringe, SGE	HP Autosampler Model 10R-HP-GT-0.47 10ul/26/R
24797	Syringe, SGE	HP Autosampler Model 10R-HP-Gt-0.63 10ul/23/R
24798	Syringe, SGE	HP Autosampler Model 25R-HP-0.63 25ul/23/R
24799	Syringe, SGE	HP Autosampler Model 50R-HP-0.63 50ul 23/R
24800	Syringe, SGE	HP Autosampler Model 100R-HP-0.63 100ul/23/R
24801	Syringe, SGE	HP Autosampler Model 250R-HP-0.63 250ul/23/R
24802	Syr Needles, SGE	Model N5-5 5.0ul/260/26 Pack of 5
24803	Syr Needles, SGE	Model N10-5 10ul/26 Pack of 5
24804	Syr Needles, SGE	Model N10-5/0.63 10ul/23 Pack of 5
24805	Syr Needles, SGE	Model 25-500ul 25 Pack of 5
24806	Syr Needles, SGE	25-500ul 23 Pack of 5
24807	Syr Needles, SGE	Model NM1/2.5-5 Pack of 5
24808	Syringe Needle, SGE	10ul 22 Gauge, 5-pk
24809	Syringe Needle, SGE	22S Gauge Pack of 5 25ul - 500ul
24810	Syringe, SGE	Model .5BR-PE-0.47 .5ul/26/R
24811	Syringe, SGE	Model .5BR-PE-0.63 .5ul/23/R
24812	Syringe, SGE	Model 5F-PE-0.47 5ul/26/F
24813	Syringe, SGE	Model 5F-PE-0.63 5ul/23/F
24814	Syringe, SGE	Model 5F-PE-GT-0.47 5ul/26/F
24815	Syringe, SGE	Model 5F-PE-GT-0.63 5ul/23/F
24816	Syringe, SGE	Model 50F-PE-0.63 50ul/23/F
24817	Syringe Needles, SGE	Model NP0.5B-PE-0.47 .5ul/26
24818	Syringe Needles, SGE	Model NP0.5B-PE-0.63 .5ul/23
24819	Inland 45 Rough Pump Oil	for MS Pumps, 1 Liter
24820	Syringe Needles, SGE	Model NP0.5B-HP-0.47 .5ul/26
24821	Syringe Needles, SGE	Model NP0.5B-HP-0.63 .5ul/23
24822	Syringe Needles, SGE	Model N5-HP-0.47 5ul/26 Pack of 2
24823	Syringe Needles, SGE	Model N5-HP-0.63 5ul/23 Pack of 2
24824	Syringe Needles, SGE	Model N10-HP-0.47 10ul/26 Pack of 2
24825	Syringe Needles, SGE	Model N10-HP-0.63 10ul/23 Pack of 2
24826	Syringe Needles, SGE	Model N25/500-HP-0.63 Pack of 2
24827	Ferrules, Graphite	0.4mm ID for 0.25 Columns and Smaller, for Shimadzu GCs, 10-pk
24828	Ferrules, Graphite	0.5mm ID for 0.32 Columns, for Shimadzu GCs, 10-pk
24829	Ferrules, Graphite	0.8mm ID for 0.53 Columns, for Shimadzu GCs, 10-pk
24830	Syringe Needles, SGE	Model NP0.5B-7C .5ul/23
24831	Syringe Needles, SGE	Model NP0.5B-7BV .5ul/23
24832	Syringe Needles, SGE	ModelNP0.5B-OC-7/0.47 .5ul/26
24833	Syringe Needles, SGE	Model NP1B-5C 1.0ul/23
24834	Syringe Needles, SGE	Model NP1B-7C 1.0ul/23
24835	Syringe Needles, SGE	Model NP1B-7BV 1.0ul/23
24836	Syringe Needles, SGE	Model NP1B-OC-7/0.47 1.0ul/26
24837	Syringe Needle, SGE	Model NP5B-5 5ul/23

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

**HROMalytic** +61(0)3 9762 2034
ECHnology Pty Ltd
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

(Ctrl F)
> Cat.#s

> **Current System AUD Prices**

