

24838 Syringe Needles, SGE
24839 Syringe Needles, SGE
24840 Syringe, SGE
24841 Syringe, SGE
24842 Syringe, SGE
24843 Syringe, SGE
24844 Syringe, SGE
24845 Syringe, SGE
24846 Syringe, SGE
24847 Syringe, SGE
24848 Ferrule, Vespel/Graphite
24849 Ferrule, Vespel/Graphite
24850 Syringe, SGE
24851 Syringe, SGE
24852 Syringe, SGE
24853 Syringe, SGE
24854 Syringe Needles, SGE
24855 Syringe Needles, SGE
24856 Syringe Needles, SGE
24857 Syringe Needles, SGE
24858 Syringe Needles, SGE
24859 Syringe Needles, SGE
24860 Syringe, SGE
24861 Syringe, SGE
24862 Syringe, SGE
24863 Syringe, SGE
24864 Syringe, SGE
24865 Syringe, SGE
24866 Syringe, SGE
24867 Syringe, SGE
24868 Syringe, SGE
24869 Syringe, SGE
24870 Syringe, SGE
24871 Syringe, SGE
24878 Syringe Needles, SGE
24879 Syringe Needles, SGE
24880 Syringe, SGE
24881 Syringe, SGE
24882 Syringe, SGE
24883 Syringe, SGE
24884 Syringe, SGE
24885 Syringe, SGE
24886 Syringe, SGE
24887 Syringe, SGE
24888 Syringe, SGE
24889 Syringe, SGE
24890 Syringe, SGE
24891 Syringe, SGE
24892 Syringe, SGE
24894 Syringe, Hamilton
24895 Syringe, Hamilton
24896 Nut for Terminal Fitting
24898 Syringe Valve
24899 Viton O-rings
24900 Miniert Valve
24901 Miniert Valve
24902 Miniert Valve
24903 Miniert Valve
24904 Miniert Valve
24905 Miniert Valve
24906 Mininert Septa
24907 Septum Insertion Tool
24912 Ferrule, Vespel/Graphite
24913 Generator, Hydrogen & Air
24913-550 Generator, Hydrogen & Air
24914 Resin Bed Cartridge

Page 108
Model NP5B-7 5.0ul/23
Model NP5B-7BV 5.0ul/23
Shimadzu Autosampler Model .5BR-S 0.17-0.63 5ul/23/R
SGE Shimadzu Autosampler Syr Model .5BR-S-0.63 .5ul/23/R
SGE Shimadzu Autosampler Syr Model 5F-S-0.47 5ul/26/F
SGE Shimadzu Autosampler Syr Model 5F-S-0.63 5ul/23/F
SGE Shimadzu Autosampler Syr Model 10R-S-0.47 10ul/26/R
SGE Shimadzu Autosampler Syr Model 10R-S-0.63 10ul/23/R
SGE Shimadzu Autosampler Syr Model 10R-S-GT-0.63 10ul/23/R
SGE Shimadzu Autosampler Syr Model 250F-S-GT-0.63 250ul/23
2-Hole, 1/16" x 0.4mm 5pk
2-Hole, 1/16" x 0.5mm 5pk
SGE Varian Autosampler Syringe Model 10r-VA80-II 10ul 25/term
SGE Varian Autosampler Syringe Model 10F-VA8X-11 10ul 25/F
SGE Varian Autosampler Syringe Model 10R-VA8X-II 10ul/25/R
Varian Autosampler Model 100R-VA8X 100ul/25/R
Model N10-VA8X00H-II 10ul/25
Model N10-VA8035-II 10ul/25 Pack of 2
ModelN10-VA8035-0.17-II 10ul Pack of 2
ModelN10-VA8100H-0.63-II 10ul 2pk
Model N10-VA800H-II(0.2) 10ul
Model N25/500-VA 100ul/25
For Rheodyne & Valco Model 10F-LC 10ul/22s/3
For Rheodyne & Valco Model 25F-LC 25ul/22s/3
For Rheodyne & Valco Model 50F-LC 50ul/22s/3
For Rheodyne & Valco Model 100F-LC 100ul/22s/3
For Rheodyne & Valco Model 250F-LC 250ul/22s/3
For Rheodyne & Valco Model 500F-LC 500ul/22s/3
GasTght Syr/Rheodyne&Valco Model 10R-GT-LC 10ul/22/3
GasTght Syr/Rheodyne&Valco Model 25R-GT-LC 25ul/22s/3
GasTght Syr/Rheodyne&Valco Model 50R-GT-LC 50ul/22s/3
GasTght Syr/Rheodyne&Valco Model 100R-GT-LC 100ul/22/3
GasTght Syr/Rheodyne&Valco Model 250R-GT-LC 250ul/22s/3
GasTght Syr/Rheodyne&Valco Model 500R-GT-LC 500ul/22s/3
Model N10-WLC 25guage 5PK
Model n25/500-WLC 25gauge 5PK
Model 50R-V-GT 50ul/23
Model 100R-V-GT 100ul/23
Model 250R-V-GT 250ul/23
Model 500R-V-GT 500ul/23
Model 1MDR-V-GT 1ml/23
Model 2.5MDR-V-GT 2.5ml/23
Model 1MR-VLL-GT 1ml/Leur Lock
Model 2.5MDR-VLL-GT 2.5ml/LL
Model 5MDR-VLLMA-GT 5ml/LeurLo
Model 10MDR-VLLMA-GT 10ml/LL
Model 25MDR-VLLMA-GT 25ml/LL
Model 50MR-VLLMA-GT 50ml/LL
Model 100MR-VLLMA-GT 100ml
Model 1701 10ul N/F GT 23s gauge
Model 7000-5 0.5ul R
for Thermo Trace GCs, 2pk
Model V fits SGE RN 25ul-2.5mL
for Shimadzu 17A, 2010 and 2014 GCs, 10pk
For 13mm Screw Cap Vials Pack of 12
For 15mm Screw Cap Vials Pack of 12
For 18mm Screw Cap Vials Pack of 12
For 20mm Screw Cap Vials Pack of 12
For 24mm Screw Cap Vials Pack of 12
For 20mm Crimp Top Vials Pack of 12
Replacement Septa Pack of 50
1/16" x 1.0mm ID Pack of 10
FID-2500 Gas Station Domestic Orders Drop Ship From Manufacturer Ships Via I
FID-2500 Gas Station United Kingdom-230 VAC 50/50 Hz 12 month warranty fron
for FID-1000 and FID-2500 Hydrogen Generators Parker/Balston

(Ctrl F)
> Cat.#s

> Current System AUD Prices



Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

**HROMalytic** +61(0)3 9762 2034
ECHnology Pty Ltd
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au