



25221	Secure-Fit Fitting	PEEK Double Ended 10cm x 0.005"
25222	Secure-Fit Fitting	PEEK Double Ended 20cm x 0.005"
25223	Secure-Fit Fitting	PEEK Double Ended 30cm x 0.005"
25224	Secure-Fit Fitting	PEEK Double Ended 10cm x 0.007"
25225	Secure-Fit Fitting	PEEK Double Ended 20cm x 0.007"
25226	Secure-Fit Fitting	PEEK Double Ended 30cm x 0.007"
25227	Secure-Fit Fitting	PEEK Double Ended 10cm x 0.010"
25228	Secure-Fit Fitting	PEEK Double Ended 20cm x 0.010"
25229	Secure-Fit Fitting	PEEK Double Ended 30cm x 0.010"
25230	Secure-Fit Fitting	PEEK Single End 6cm x 0.005"
25231	Secure-Fit Fitting	PEEK Single End 10cm x 0.005"
25232	Secure-Fit Fitting	PEEK Single End 20cm x 0.005"
25233	Secure-Fit Fitting	PEEK Single End 30cm x 0.005"
25234	Secure-Fit Fitting	PEEK Single End 6cm x 0.007"
25235	Secure-Fit Fitting	PEEK Single End 10cm x 0.007"
25236	Secure-Fit Fitting	PEEK Single End 20cm x 0.007"
25237	Secure-Fit Fitting	PEEK Single End 30cm x 0.007"
25238	MiniPulse Pulse Damper	160uL Dead Volume 316 Stainless Steel
25239	MiniPulse Pulse Damper	160uL Dead Volume, PEEK
25240	Capillary Tubing, SS	1/16"OD x 0.005 ID x 5cm L, 3pk
25241	Capillary Tubing, SS	1/16"OD x 0.005 ID x 10cm L, 3pk
25242	Capillary Tubing, SS	1/16"OD x 0.005 ID x 20cm L, 3pk
25243	Capillary Tubing, SS	1/16"OD x 0.005 ID x 30cm L, 3pk
25244	Capillary Tubing, SS	1/16"OD x 0.007 ID x 5cm L, 3pk
25245	Capillary Tubing, SS	1/16"OD x 0.007 ID x 10cm L, 3pk
25246	Capillary Tubing, SS	1/16" OD x 0.007 ID x 20cm L, 3pk
25247	Capillary Tubing, SS	1/16" OD x 0.007 ID x 30cm L, 3pk
25248	Capillary Tubing, SS	1/16" OD x 0.01 ID x 5cm L, 3pk
25249	Capillary Tubing, SS	1/16"OD x 0.01 ID x 10cm L, 3pk
25250	Capillary Tubing, SS	1/16" OD x 0.01 ID x 20cm L, 3pk
25251	Capillary Tubing, SS	1/16" OD x 0.01 ID x 30cm L, 3pk
25252	Capillary Tubing, SS	1/16"OD x 0.02 ID x 5cm L, 3pk
25253	Capillary Tubing, SS	1/16"OD x 0.02 ID x 10cm L, 3pk
25254	Capillary Tubing, SS	1/16"OD x 0.02 ID x 20cm L, 3pk
25255	Capillary Tubing, SS	1/16"OD x 0.02 ID x 30cm L, 3pk
25256	PEEK Washer	for ASE 200 Extraction Unit 12pk, Replaces Dionex Part #049454
25257	PEEK Washer	for ASE 200 Extraction Unit 48pk, Replaces Dionex Part #049454
25258	Needle Seat	Agilent 1050
25259	Preventative Maint. Kit	for Agilent 1050 Autosampler Similar to Agilent #01078-68721
25260	Deuterium Lamp	Agilent 1050 VW/DA/MW, 1090 DA, 1040/1040M
25261	Deuterium Lamp	1100 DAD G1315A, G1365A, Agilent 1100/1200
25262	Deuterium Lamp	Agilent 1100/1200 VWD
25263	Deuterium Lamp	Agilent 8453
25264	Lamp	G1321 Fluorescence Detector Flash Lamp
25265	Needle Seat Assembly	Agilent 1100
25266	Sieves for Outlet Valve	Agilent 1050/1100/1200, 10pk
25267	Outlet Ball Valve	Binary Pump, Agilent 1100/1200
25268	Seal Wash Kit	Binary Pump, Agilent 1100 4 Seals, 4 Gaskets
25269	Seal Wash Kit	2 Seals, 2 Gaskets Agilent 1050/1100
25270	Pump Maintenance Kit	Agilent 1050/1100
25271	Maintenance Kit	Autosampler Preventative Maint. Kit, Agilent 1100
25272	Rotor Seal	Agilent 1050, Not for use with 7125
25273	Sapphire Piston	Agilent 1050/1100/1200
25275	Rotor Seal	Agilent 1100
25276	Outlet Ball Valve	Agilent 1050/1100/1200
25277	Wash Seal	Agilent 1050/1100/1200
25278	Needle Assembly	Agilent 1100
25279	7010 Vespel Rotor Seal	
25280	7125 Vespel Rotor Seal	
25281	Isolation Seal, 7010	Rheodyne Style
25282	Outlet Check Valve	LC-600, LC-9A, LC-10AD, 10AT
25284	Deuterium Lamp	Shimadzu SPD-10A, 10AVP 10AVVP
25285	Plunger Seal	LC-6A
25287	Inlet Check Valve	LC-6A, LC-10AS
25288	Outlet Check Valve	LC-6A, LC-10AS
25289	Check Valve Rebuild Kit	for Outlet Check Valve 2pk,LC-6A, LC-10AS