

25640 GC Packing Material  
 25641 INDEFINITE BACKORDER  
 25642 GC Packing Material  
 25643 GC Packing Material  
 25644 GC Packing Material  
 25645 GC Packing Material  
 25646 GC Packing Material  
 25647 GC Packing Material  
 25648 GC Packing Material  
 25649 GC Packing Material  
 25650 GC Packing Material  
 25651 GC Packing Material  
 25652 GC Packing Material  
 25653 GC Packing Material  
 25654 GC Packing Material  
 25655 GC Packing Material  
 25656 GC Packing Material  
 25657 GC Packing Material  
 25660 GC Packing Material  
 25661 GC Packing Material  
 25664 GC Packing Material  
 25665 GC Packing Material  
 25668 GC Packing Material  
 25672 GC Packing Material  
 25675 GC Packing Material  
 25676 GC Packing Material  
 25678 GC Packing Material  
 25693 GC Packing Material  
 25700 GC Packing Material  
 25701 GC Packing Material  
 25959 Replacement Frits  
 25988 SPE C18-90mm Disks  
 25989 DVB SPE Tubes, 6mL, 500mg  
 25990 Cap Frits 2mm 0.5um  
 25991 Glass Vacuum Chamber  
 25992 Resprep Disk Manifold  
 26000 Endcap, SPE Tube, 12pk  
 26001 SPE Lid 1ml Tube, 12pk  
 26002 SPE Lid 3ml Tube, 12pk  
 26003 SPE Lid 6ml Tube, 12pk  
 26004 SPE Lid 15ml Tube, 12pk  
 26005 SPE Lid 25ml Tube, 12pk  
 26006 SPE Lid 70ml Tube, 12pk  
 26007 SPE Connector 1,3,6ml  
 26008 SPE Connector 25ml  
 26009 SPE Connector 60ml  
 26010 SPE empty 1ml Tubes,50pk  
 26011 SPE empty 3ml Tubes,50pk  
 26012 SPE empty 6ml Tubes,50pk  
 26013 SPE empty 15ml Tubes,50pk  
 26014 SPE sample reservoir 25ml  
 26015 SPE sample reservoir 70ml  
 26016 SPE PE Frits 1ml Tubes  
 26017 SPE PE Frits 3ml Tubes  
 26018 SPE PE Frits 6ml Tubes  
 26019 SPE PE Frits 15ml Tubes  
 26020 SPE PE Frits 25ml Tubes  
 26021 SPE PE Frits 70ml Tubes  
 26022 SPE Oil&Grease Fiber Disk  
 26023 SPE Resin-47mmdisks,pk 20  
 26024 SPE Flowfilter-47mm,pk 20  
 26029 AH SPE Tubes 6mL, 500mg  
 26030 SPE C18,1mL/100mg,100pk  
 26031 SPE C18,3mL/200mg, 50pk  
 26032 SPE EPA Method 521 & 522  
 26033 Diatomaceous Earth

Page 118  
 INDEFINITE BACKORDER One Bottle  
 Silcoport P 80/100 MINIMUM ORDER 100 GRAMS  
 Silcoport P 100END20INDEFINITE B.  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams INDEFINITE BACKORDER  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER One Bottle of 100 Grams (Pre-2009 sold by the gm w/  
 INDEFINITE BACKORDER Chromosorb WNAW 45/60 One Bottle of 100 Grams  
 INDEFINITE BACKORDER Chromosorb WNAW 60/80 One Bottle of 100 Grams  
 INDEFINITE BACKORDER Chromosorb WAW 45/60 One Bottle of 100 Grams (F  
 INDEFINITE BACKORDER Chromosorb WAW 60/80 One Bottle of 100 Grams (F  
 INDEFINITE BACKORDER Chromosorb WAW/DMDCS 45/60 One Bottle of 100  
 INDEFINITE BACKORDER Chromosorb WAW/DMDCS 60/80 One Bottle of 100  
 INDEFINITE BACKORDER Chromosorb W-HP 60/80 One Bottle of 100 Grams (F  
 Silcoport W 80/100 MINIMUM ORDER 100 GRAMS PACKING MATERIAL ONLY  
 Chromosorb T 30/60 MINIMUM ORDER 100 GRAMS PACKING MATERIAL ONL  
 Chromosorb T 40/60 One Bottle of 100 Grams (Pre-2009 sold by the gm w/100gr  
 High Carbon n-octane on ResSil C 80/100, 10gm Minimum Order  
 Chromosorb T 40/60 One Bottle of 100 Grams (Pre-2009 sold by the gm w/100gr  
 Tenax-TA 20/35 Min Order 10 grams  
 Tenax-TA 35/60 Min Order 10 grams  
 For ASE 200 Extraction Unit, 100-pk  
 Pack of 12  
 Teflon Frits, 50pk  
 Replacement for Trident Guard System, Pack of 5  
 w/ gauge & valve, Replacement for Resprep 6 & 12-port manifolds  
 6-port, Lid & Collection Plates, fits Resprep 12-port glass chamber (not included)  
 Universal female luer plug  
 12pk  
 12pk  
 12pk  
 12pk  
 100pk,6mm polyethylene  
 100pk,9mm polyethylene  
 100pk,1.2cm polyethylene.513"+-.005OD x .062"+-.01  
 100pk,1.6cm polyethylene  
 100pk,2.0cm polyethylene  
 100pk,2.6cm polyethylene  
 INDEFINITE BACKORDER 47mm 20pack  
 styrene/DVB resin  
 for Phenoxy Acid Herbicides 30pk, EPA Method 8321  
 6mL/2g Activated Charcoal, 30pk  
 1kg, Replaces Dionex PN 062819

