



26034	SPE C18,20mL/5g, 20pk	
26035	SPE C18,60mL/10g, 16pk	
26036	SPE C8,1mL/100mg,100pk	
26037	SPE C8,3mL/200mg, 50pk	
26038	SPE C8,3mL/500mg, 50pk	
26039	SPE C8,6mL/500mg, 30pk	
26040	SPE C8,6mL/1g, 30pk	
26041	SPE C2,1mL/100mg, 100pk	
26042	SPE C2,3mL/200mg, 50pk	
26043	SPE Cyclohexyl,6mL/1g	30pk
26044	SPE DrugPrepl,3mL/200mg	50pk,cation exchanger
26045	SPE DrugPrepl,10mL/200mg	50pk,cation exchanger
26046	SPE DrugPreplI,3mL/200mg	50pk,anion exchanger
26047	SPE DrugPreplI,10mL/200mg	50pk,anion exchanger
26048	SPE CN,3mL/200mg, 50pk	
26049	SPE CN,3mL/500mg, 50pk	
26050	SPE NH2,1mL/100mg, 100pk	
26051	SPE NH2,3mL/200mg, 50pk	
26052	SPE NH2,3mL/500mg, 50pk	
26053	SPE NH2,6mL/500mg, 30pk	
26054	SPE SAX,1mL/100mg, 100pk	
26055	SPE SAX,3mL/500mg, 50pk	
26056	SPE SCX,1mL/100mg, 100pk	
26057	SPE SCX,3mL/500mg, 50pk	
26058	SPE BenzSCX,3mL/200mg	50pk
26059	SPE BenzSCX,6mL/500mg	30pk
26060	SPE BenzSCX,6mL/1g	30pk
26061	SPE COOH,1mL/100mg, 100pk	
26062	SPE COOH,3mL/500mg, 50pk	
26063	SPE EPA 548.1,6mL, 30pk	BioRex resin, slurry pack
26065	SPE EPH Fractionation	20mL, 5g Si, 20-pk
26077	Vac.manifold,12port	EA, complete
26078	Vac.manifold,12port	EA,replacement lid w/valves
26079	Vac.manifold,12port	EA,replacement collection rack
26080	Vac.manifold,24port	EA, Complete
26081	Vac.manifold,24port	EA,replacement lid w/valves
26082	Vac.manifold,24port	EA,replacement collection rack
26083	Vac.manifold replacement	2pk, sample valves
26084	Vac.manifold replacement	12pk, Tef. sample guides
26085	SPE Cartridge	Florisil Packing Glass Tube with Teflon Frits 6ml/1g, Pack of 30.
26086	SPE Cartridge	Florisil Packing Glass Tube with Teflon Frits 6ml/500mg, Pack of 30.
26087	SPE Cartridge	Carboprep 200 Packing 6ml, 500mg, 30pk 200 m2/gm surface
26088	SPE Cartridge	Carboprep 200 Packing 3ml, 250mg 50pk 200 m2/gm surface
26091	SPE Cartridge	CarboPrep 90 Packing 3ml, 250mg 50pk 90 m2/gm surface
26092	SPE Cartridge	CarboPrep 90 Packing 6ml, 500mg 30pk 90 m2/gm surface
26093	SPE Tubes	Resprep RDX 6ml/500mg 30 pack
26096	Replacement Caps	for ASE 200 Extraction Unit 2-pk
26098	22mL Extraction Cell Body	for ASE 200 Extraction Unit
26100	Replacement Frits	for ASE 200 Extraction Unit 10pk
26110	1ml Extraction Cell Body	for ASE 200 Extraction Unit
26112	5ml Extraction Cell Body	for ASE 200 Extraction Unit
26114	11ml Extraction Cell Body	for ASE 200 Extraction Unit
26116	33ml Extraction Cell Body	for ASE 200 Extraction Unit
26118	Cellulose Filters	for ASE 200, Pack of 100 Replaces Dionex PN 049458
26119	Glass Fiber Filters	for ASE 200, Pack of 100 Replaces Dionex PN 047017
26120	PEEK Washer	for ASE 200 Extraction Unit, 250pk Replaces Dionex #049454
26121	60ml Collection Vial	Pre-Cleaned Clear Glass for ASE System, 72pk, Replaces Dionex #048784
26122	60ml Collection Vial	Precleaned Amber Glass for ASE System, 72pk. Replaces Dionex #048781
26123	QuEChERS SPE	Q-sep Q253 2mL Micro-Centrifuge Tube, 150mg MgSO ₄ , 50mg PSA, 50mg Carb
26124	QuEChERS SPE	Q-sep Q250 2mL Centrifuge Tube,150mg MgSO ₄ , 50mg PSA, 100pk Packaged i
26125	QuEChERS SPE	Q-sep Q251 Micro-Centrifuge Tube 2mL, 150mg MgSO ₄ , 50mg PSA, 50mg C18,
26126	QuEChERS SPE	Q-sep Q374 15mL Centrifuge Tube, 900mg MgSO ₄ , 300mg PSA, 150mg Carbon
26127	Combo SPE Tube	6ml SPE Tube Packed w/200mg Graphitized Carbon and 400mg PSA Teflon Frits
26128	Combo SPE Tube	6ml SPE Tube Packed w/250mg Graphitized Carbon and 500mg PSA Teflon Frits
26129	Combo SPE Tube	6ml SPE Tube Packed w/500mg Graphitized Carbon and 500mg PSA Teflon Frits
26130	SPE Manifold Part	Valve Attachment, 48pk