

26131	SPE Manifold Part	Needle Attachment, 48pk
26132	ASE 300 Tube, 100ml	Replacement Extraction Tube
26134	ASE 300/100 Snap Ring	for Caps, 12pk
26135	Florisil PR	500 grams, 60/100 mesh
26136	Activated Copper	Granulated 30 mesh, 1kg
26137	Ottawa Sand	5kg, 20/30 mesh
26138	Sparger, Fritted	10mL, for Tekmar 2000, 3000, and 3100
26139	Sparger, Non-Fritted	5mL, for Tekmar 2000, 3000, and 3100
26140	Sparger, Non-Fritted	10mL, for Tekmar 2000, 3000, and 3100
26141	Sparger, Non-Fritted	25mL, for Tekmar 2000, 3000, and 3100
26142	Syringe Filter, 13mm	0.22um PTFE 100pk Purple Luer Lock Inlet
26142	Syringe Filter, 13mm	0.22um PTFE 100pk Purple Luer Lock Inlet
26142	Syringe Filter, 13mm	0.22um PTFE 100pk Purple Luer Lock Inlet
26143	Syringe Filter, 13mm	0.45um PTFE 100pk Orange Luer Lock Inlet
26143	Syringe Filter, 13mm	0.45um PTFE 100pk Orange Luer Lock Inlet
26143	Syringe Filter, 13mm	0.45um PTFE 100pk Orange Luer Lock Inlet
26144	Syringe Filter, 25mm	0.22um PTFE 100pk Purple Luer Lock Inlet
26144	Syringe Filter, 25mm	0.22um PTFE 100pk Purple Luer Lock Inlet
26144	Syringe Filter, 25mm	0.22um PTFE 100pk Purple Luer Lock Inlet
26145	Syringe Filter, 25mm	0.45um PTFE 100pk Orange Luer Lock Inlet
26145	Syringe Filter, 25mm	0.45um PTFE 100pk Orange Luer Lock Inlet
26145	Syringe Filter, 25mm	0.45um PTFE 100pk Orange Luer Lock Inlet
26146	Syringe Filter, 13mm	0.22um Nylon 100pk Luer Lock Inlet, Yellow
26146	Syringe Filter, 13mm	0.22um Nylon 100pk Luer Lock Inlet, Yellow
26146	Syringe Filter, 13mm	0.22um Nylon 100pk Luer Lock Inlet, Yellow
26147	Syringe Filter, 13mm	0.45um Nylon 100pk Luer Lock Inlet, Pink
26147	Syringe Filter, 13mm	0.45um Nylon 100pk Luer Lock Inlet, Pink
26147	Syringe Filter, 13mm	0.45um Nylon 100pk Luer Lock Inlet, Pink
26148	Syringe Filter, 25mm	0.22um Nylon 100pk Yellow Luer Lock Inlet
26148	Syringe Filter, 25mm	0.22um Nylon 100pk Yellow Luer Lock Inlet
26148	Syringe Filter, 25mm	0.22um Nylon 100pk Yellow Luer Lock Inlet
26149	Syringe Filter, 25mm	0.45um Nylon 100pk Luer Lock Inlet, Pink
26149	Syringe Filter, 25mm	0.45um Nylon 100pk Luer Lock Inlet, Pink
26149	Syringe Filter, 25mm	0.45um Nylon 100pk Luer Lock Inlet, Pink
26150	Syringe Filter, 13mm	0.22um PVDF 100pk Brown Luer Lock Inlet
26150	Syringe Filter, 13mm	0.22um PVDF 100pk Brown Luer Lock Inlet
26150	Syringe Filter, 13mm	0.22um PVDF 100pk Brown Luer Lock Inlet
26151	Syringe Filter, 13mm	0.45um PVDF 100pk Red Luer Lock Inlet
26151	Syringe Filter, 13mm	0.45um PVDF 100pk Red Luer Lock Inlet
26151	Syringe Filter, 13mm	0.45um PVDF 100pk Red Luer Lock Inlet
26152	Syringe Filter, 25mm	0.22um PVDF 100pk Brown Luer Lock Inlet
26152	Syringe Filter, 25mm	0.22um PVDF 100pk Brown Luer Lock Inlet
26152	Syringe Filter, 25mm	0.22um PVDF 100pk Brown Luer Lock Inlet
26153	Syringe Filter, 25mm	0.45um PVDF 100pk Red Luer Lock Inlet
26153	Syringe Filter, 25mm	0.45um PVDF 100pk Red Luer Lock Inlet
26153	Syringe Filter, 25mm	0.45um PVDF 100pk Red Luer Lock Inlet
26155	Syringe Filter, 13mm	0.45um Cellulose Acetate 100pk blue, luer lock
26155	Syringe Filter, 13mm	0.45um Cellulose Acetate 100pk blue, luer lock
26155	Syringe Filter, 13mm	0.45um Cellulose Acetate 100pk blue, luer lock
26156	Syringe Filter, 13mm	0.22um Cellulose Acetate 100pk Green, luer lock
26156	Syringe Filter, 13mm	0.22um Cellulose Acetate 100pk Green, luer lock
26156	Syringe Filter, 13mm	0.22um Cellulose Acetate 100pk Green, luer lock
26157	Syringe Filter, 25mm	0.45um Cellulose Acetate 100pk Blue, luer lock
26157	Syringe Filter, 25mm	0.45um Cellulose Acetate 100pk Blue, luer lock
26157	Syringe Filter, 25mm	0.45um Cellulose Acetate 100pk Blue, luer lock
26158	Syringe Filter, 25mm	0.22um Cellulose Acetate 100pk Green, luer lock
26158	Syringe Filter, 25mm	0.22um Cellulose Acetate 100pk Green, luer lock
26158	Syringe Filter, 25mm	0.22um Cellulose Acetate 100pk Green, luer lock
26159	Pressure Relief Valve	for ASE 100, 200 and 300 Systems
26160	Static Valve	for ASE 100, 200 and 300 Systems
26161	Autoseal Tip Assembly	for ASE 200 Systems, Siltek Treated
26162	Autoseal Tip Assembly	for ASE 200 Systems, Stainless Steel
26166	ASE 200 Cap Insert	for use with Cat. #26096 only 2pk
26167	ASE 300 Cap Insert	for use with Cat. #26170 only 2pk
26168	Cellulose Filters	for ASE 300, 100pk Replaces Dionex #056780
26169	ASE Funnel	For ASE 300 Sample Cells

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

		Australian Distributors Importers & Manufacturers www.chromtech.net.au
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA		

(Ctrl F)
> Cat.#s

> Current System AUD Prices

