

26170 ASE Caps
 26172 ASE 300 Tube, 10ml
 26174 ASE Frits
 26176 ASE 300 Tube, 34ml
 26178 ASE 300 Tube, 66ml
 26180 ASE Tools
 26181 ASE Tools
 26182 ASE Tools
 26183 ASE Tools
 26184 ASE 200 Snap Ring
 26185 ASE Pliers
 26187 O-rings, PTFE
 26188 O-rings, Viton
 26189 Filters, Glass Fiber
 26190 Filters, Cellulose
 26191 Collection Bottle
 26193 Combo SPE Tube
 26194 Combo SPE Tube
 26195 SPE Tube
 26204 Sodium Sulfate
 26205 SPE Florisil, 6mL/1000mg
 26206 Closed End SPE Cartridges
 26207 Closed End SPE Cartridges
 26208 SPE GHB Tubes
 26209 1/4" Sample Inlet
 26210 1/4" Sample Inlet
 26211 1/4" Sample Inlet
 26212 1/4" Sample Inlet
 26215 QuEChERS SPE
 26216 QuEChERS SPE
 26217 QuEChERS SPE
 26218 QuEChERS SPE
 26219 QuEChERS SPE
 26220 QuEChERS SPE
 26221 QuEChERS SPE
 26222 QuEChERS SPE
 26223 QuEChERS SPE
 26224 QuEChERS SPE
 26225 QuEChERS SPE
 26226 QuEChERS SPE
 26228 SPE Tubes
 26229 PEEK Washer
 26230 Q-sep 3000 Centrifuge
 26231 Q-sep 3000 Centrifuge
 26232 50mL Tube Carrier
 26233 4-place Tube Carrier
 26234 2mL Tube Adapters
 26235 QuEChERS SPE
 26236 QuEChERS SPE
 26237 QuEChERS SPE
 26238 QuEChERS SPE
 26239 QuEChERS SPE
 26246 Autoseal Tip Assembly
 26247 Autoseal Tip Assembly
 26248 Tubing Assembly
 26249 50mL Conical Tube Insert
 26250 Sample Delivery System
 26251 Tubing Assembly
 26252 Mini Passive Air Sampling
 26253 Mini Passive Air Sampling
 26254 Mini Passive Air Sampling
 26255 Mini Passive Air Sampling
 26256 Mini Passive Air Sampling
 26257 Mini Passive Air Sampling
 26258 Mini Passive Air Sampling
 26259 Mini Passive Air Sampling

Page 121
 ASE 300 Extraction Replacement Cap
 Replacement Extraction Tube
 ASE 300/100 Replacement Frits, 6pk
 Replacement Extraction Tube
 Replacement Extraction Tube
 Funnel for ASE 200
 2-in-1 Filter/O-Ring Insertion Tool for ASE 100/200/300 (includes Resprep Tool H:
 Resprep Tool Handle for ASE 100/200/300
 Filter Insertion Attachments for ASE 100/200/300 (1mL, 5mL, 11mL, 33mL)
 for caps, 10pk
 ASE Internal Retaining Ring Pliers, for use with ASE 200 and 300 Sample Tubes
 for ASE 200 and 300 Caps, 100pk, Replaces Dionex #049457
 for ASE 200 and 300 Caps, 50pk Replaces Dionex #056325
 for ASE 300, 100pk, Replaces Dionex #056781
 for ASE 200, 1000pk, Replaces Dionex #049458
 250mL, Pre-cleaned Clear Glass for ASE 300, 12pk, Replaces Dionex #056284
 6mL/500mg Carboprep 90/500 Aminopropyl Polyethylene Frits 30pk
 6mL/500mg Carboprep 90/500mg PSA Polyethylene Frits 30pk
 90/500mg PSA Polyethylene Frits 30pk
 5kg
 100pk
 Carboprep 200, 500mg, 30pk
 Activated Sodium Sulfate 2g, 50pk
 200mg, 10ml, 50pk
 w/Nut and Ferrules, 6" Length
 w/Nut and Ferrules, 6" Length, Siltek Treated
 w/Nut and Ferrules, 1.5" Length
 w/Nut and Ferrules, 1.5" Length Siltek Treated
 Q-sep Q210 Micro-Centrifuge Tube 2mL 150mg, MgSO₄, 25mg PSA, 100-pk Pack
 Q-sep Q211 2mL Micro-Centrifuge Tube, 150mg MgSO₄, 25mg PSA, 25mg C18,
 Q-sep Q212 Micro-Centrifuge Tube, 2mL 150mg MgSO₄, 25mg PSA, 2.5mg Car
 Q-sep Q213 Micro-Centrifuge Tube, 2mL, 150mg MgSO₄, 25mg PSA, 7.5mg Car
 Q-sep Q252 Micro-Centrifuge Tube 2mL, 150mg MgSO₄, 50mg PSA, 50mg C18,
 Q-sep Q350 15mL Centrifuge Tube, 1200mg MgSO₄, 400mg PSA, 50pk, Packag
 Q-sep Q351, 15mL Centrifuge Tube, 1200mg MgSO₄, 400mg PSA, 400mg C18, !
 Q-sep Q352 15mL Centrifuge Tube, 1200mg MgSO₄, 400mg PSA, 400mg C18, 4
 Q-sep Q370 15mL Centrifuge Tube, 900mg MgSO₄, 150mg PSA, 50pk
 Q-sep Q371 15mL Centrifuge Tube, 900mg MgSO₄, 150mg PSA, 15mg Carbon, !
 Q-sep Q372 15mL Centrifuge Tube 900mg MgSO₄, 150mg PSA, 45mg Carbon, 5
 Q-sep Q373 15mL Centrifuge Tube, 900mg MgSO₄, 150mg PSA, 150mg C18, 50
 15mL/2g Florisil, 100-pk
 for ASE 200 Extraction Unit, 1000-pk, Replaces Dionex # 049454
 110V
 220V
 for the Q-sep 3000 Centrifuge, 2-pk
 for the Q-sep 3000 Centrifuge, 2-pk
 for the Q-sep 3000 Centrifuge, 4-pk
 Q-sep Q110 Packet & 50mL Tube, 4gMgSO₄, 1gNaCl, 1gTSCD, 0.5gDHS, Box of !
 Q-sep Q110 Packet 4gMgSO₄, 1gNaCl, 1gTSCD, 0.5gDHS Box of 50 Packets
 Q-sep Q150 Packet & 50ml Tube 6g MgSO₄, 1.5g NaOAc Box of 50 Packets and
 Q-sep Q150 Packet 6g MgSO₄, 1.5g NaOAc Box of 50 Packets
 Q-sep Empty 50ml Tube, 50-pk
 For Dionex ASE 300 Systems, Stainless Steel
 For Dionex ASE 300 Systems, Siltek treated
 For Dionex ASE 300 Systems
 for the Q-sep 3000 Centrifuge, 6-pk
 Compatible with Resprep Tubes 1ml, 3ml, 6ml, 15ml and Diskcover-47, 6pk
 for Dionex ASE 200 Systems
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Non-Siltek, 0.000
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Siltek, 0.0008" Oi
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Non-Siltek, 0.001
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Siltek, 0.0012" Oi
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Non-Siltek, 0.001
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Siltek, 0.0016" Oi
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Non-Siltek, 0.002
 CONFIRM DELIVERY WITH OPERATIONS > 25 Complete Kit, Siltek, 0.0020" Oi

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
 Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

