

26260	Collection Bottle
26261	Tool
26262	Tool
26263	Grab Sampling Kit
26264	Grab Sampling Kit
26265	Grab Sampling Kit
26266	Grab Sampling Kit
26267	Grab Sampling Kit
26268	Grab Sampling Kit
26269	Grab Sampling Kit
26270	Grab Sampling Kit
26271	Replacement Fitting
26272	Replacement Fitting
26273	Replacement Fitting
26274	Replacement Fitting
26275	Replacement Fitting
26276	Replacement Fitting
26277	Replacement Fitting
26278	Replacement Fitting
26279	10 Micron Frit
26280	Siltek Grab Sampling Kit
26281	Siltek Grab Sampling Kit
26282	Siltek Grab Sampling Kit
26283	Siltek Grab Sampling Kit
26284	Siltek Grab Sampling Kit
26285	Siltek Grab Sampling Kit
26286	Siltek Grab Sampling Kit
26287	Siltek Grab Sampling Kit
26288	Replacement Fitting
26289	Replacement Fitting
26290	Replacement Fitting
26291	Replacement Fitting
26292	Replacement Fitting
26293	Replacement Fitting
26294	Replacement Fitting
26295	Replacement Fitting
26296	Siltek 10 micron Frit
26297	Sample Cylinder Valve
26298	Sample Cylinder Valve
26299	Sample Cylinder Valve
26300	Sample Cylinder Valve
26301	USE CTO.MISC TO PROCESS
26301-1	Sample Cylinder Valve
26301-2	Sample Cylinder Valve
26302	Sample Cylinder Valve
26303	Sample Cylinder Valve
26304	Sample Cylinder Valve
26305	Sample Cylinder Valve
26306	Sample Cylinder Valve UHP
26307	Sample Cylinder Valve UHP
26308	Sample Cylinder Valve UHP
26309	Sample Cylinder Valve UHP
26310	Sample Cylinder Valve UHP
26311	Sample Cylinder Valve UHP
26312	USE CTO.MISC TO PROCESS
26312-1	Sample Cylinder Valve
26313	USE CTO.MISC TO PROCESS
26313-1	Sample Cylinder Valve
26313-2	Sample Cylinder Valve
26314	Sample Cylinder Valve UHP
26315	Sample Cylinder Valve UHP
26316	Sample Cylinder Valve UHP
26317	Sample Cylinder Valve UHP
26318	Rupture Disc Tee
26319	Rupture Disc Tee
26320	Rupture Disc Replacement

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd
 Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
 Importers & Manufacturers
www.chromtech.net.au

250mL, Precleaned Amber Glass for ,

The Claw

The Claw Holder Kit, Includes The Claw and Frit.

Includes 6" sample inlet, fitting w/ 0.0018" orifice, 10 micron frit, Stainless Steel

Includes 6" sample inlet, fitting w/ 0.0020" orifice, 10 micron frit, Stainless Steel

Includes 6" sample inlet, fitting w/ 0.0030" orifice, 10 micron frit, Stainless Steel

Includes 6" sample inlet, fitting w/ 0.0040" orifice, 10 micron frit, Stainless Steel

Includes 6" sample inlet, fitting w/ 0.0055" orifice, 10 micron frit, Stainless Steel

Includes 6" sample inlet, fitting w/ 0.0080" orifice, 10 micron frit, Stainless Steel

Includes 6" sample inlet, fitting w/ 0.0090" orifice, 10 micron frit, Stainless Steel

Includes 6" sample inlet, fitting w/ 0.0130" orifice, 10 micron frit, Stainless Steel

w/ 0.0018" orifice, for Grab Sampling Kit, SS

w/ 0.0020" orifice, for Grab Sampling Kit, SS

w/ 0.0030" orifice, for Grab Sampling Kit, SS

w/ 0.0040" orifice, for Grab Sampling Kit, SS

w/ 0.0055" orifice, for Grab Sampling Kit, SS

w/ 0.0080" orifice, for Grab Sampling Kit, SS

w/ 0.0090" orifice, for Grab Sampling Kit, SS

w/ 0.0130" orifice, for Grab Sampling Kit, SS

Replacement for Grab Sampling Kit, Stainless Steel, 3-pk

Includes 6" sample inlet, fitting w/ 0.0018" orifice, 10 micron frit, Siltek Treated

Includes 6" sample inlet, fitting w/0.0020" orifice, 10 micron frit, Siltek Treated

Includes 6" sample inlet, fitting w/ 0.0030" orifice, 10 micron frit, Siltek Treated

Includes 6" sample inlet, fitting w/ 0.0040" orifice, 10 micron frit, Siltek Treated

Includes 6" sample inlet, fitting w/ 0.0055" orifice, 10 micron frit, stainless steel

Includes 6" sample inlet, fitting w/ 0.0080" orifice, 10 micron frit, stainless steel

Includes 6" sample inlet, fitting w/ 0.0090" orifice, 10 micron frit, stainless steel

Includes 6" sample inlet, fitting w/ 0.0130" orifice, 10 micron frit, stainless steel

w/ 0.0018" orifice, for Grab Sampling Kit, Siltek treated

w/ 0.0020" orifice, for Grab Sampling Kit, Siltek treated

w/ 0.0030" orifice, for Grab Sampling Kit, Siltek treated

w/ 0.0040" orifice, for Grab Sampling Kit, Siltek treated

w/ 0.0055" orifice, for Grab Sampling Kit, Siltek treated

w/ 0.0080" orifice, for Grab Sampling Kit, Siltek treated

w/ 0.0090" orifice, for Grab Sampling Kit, Siltek treated

w/ 0.0130" orifice, for Grab Sampling Kit, Siltek treated

Replacement for Grab Sampling Kit, Siltek Treated, 3-pk

1/4" Male NPT x 1/4" Male NPT 3500psig DOT press rating Stainless Steel

1/4" Male NPT x 1/4" Female NPT 3500psig DOT press rating, Stainless Steel

1/4" Male NPT x 1/4" Female NPT 3500psig DOT press rating, Sulfinert

1/4" Male NPT x 1/4" Male Compression, 3500psig DOT press rating, SS

USE CTO.MISC TO PROCESS ORDER Sample Cylinder Valve 1/4" Male NPT x

SS, 1/4"NPTx1/4"NPT w/Dip Tube 1-1/2 inches

SS, 1/4"NPTx1/4"NPT w/Dip Tube 3-3/8 inches

1/4" Male NPT x 1/4" Male NPT w/1800psi Rupture Disc, 3500psig DOT press rat

1/4" Male NPT x 1/4" Male NPT w/1800psi Rupture Disc, 3500psig DOT press rat

1/4" Male NPT x 1/4" Female NPT w/1800psi Rupture Disc, 3500psig DOT press

1/4" Male NPT x 1/4" Female NPT w/1800psi Rupture Disc, 3500psig DOT press

1/4" Male NPT x 1/4" Male NPT 5000psig DOT press rating, Stainless Steel

1/4" Male NPT x 1/4" Male NPT 5000psig DOT press rating, Sulfinert

1/4" Male NPT x 1/4" Female NPT 5000psig DOT press rating, Stainless Steel

1/4" Male NPT x 1/4" Female NPT 5000psig DOT press rating, Sulfinert

1/4" Male NPT x 1/4" Male Compression, 5000psig DOT press rating, SS

1/4" Male NPT x 1/4" Male Compression, 5000psig DOT press rating, Sulfinert

USE CTO.MISC TO PROCESS ORDER Sample Cylinder Valve UHP 1/4" Male N

SS, 1/4"NPTx1/4"NPT w/Dip Tube 1-1/2 inches

USE CTO.MISC TO PROCESS ORDER Sample Cylinder Valve UHP 1/4" Male N

Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 1-1/2 inches

Sulfinert, 1/4"NPTx1/4"NPT w/Dip Tube 5-1/4 inches

1/4" Male NPT x 1/4" Male NPT w/2850psi Rupture Disc, 5000psig DOT press rat

1/4" Male NPT x 1/4" Male NPT w/2850psi Rupture Disc, 5000psig DOT press rat

1/4" Male NPT x 1/4" Female NPT w/2850psi Rupture Disc, 5000psig DOT press

1/4" Male NPT x 1/4" Female NPT w/2850psi Rupture Disc, 5000psig DOT press

1/4" Male NPT x 1/4" Female NPT 1800psig DOT press rating Stainless Steel

1/4" Male NPT x 1/4" Female NPT 1800psig DOT press rating Sulfinert

Fits Rupture Disc Tee 1800psig DOT press rating Stainless Steel

