

910550210	Ultra Phenyl	10 x 4.0mm 3-pk
910550212	Ultra Phenyl	10 x 2.1mm 3-pk
910550220	Ultra Phenyl	20 x 4.0mm 2-pk
910550222	Ultra Phenyl	20 x 2.1mm 2-pk
9105510	Ultra Phenyl 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9105511	Ultra Phenyl 5um	100 x 1.0mm
9105512	Ultra Phenyl 5um	100 x 2.1mm
9105512-700	Ultra Phenyl 5um	100 x 2.1mm with Trident
9105513	Ultra Phenyl 5um	100 x 3.2mm
9105513-700	Ultra Phenyl 5um	100 x 3.2mm with Trident
9105515	Ultra Phenyl 5um	100 x 4.6mm
9105515-700	Ultra Phenyl 5um	100 x 4.6mm with Trident
9105517	Ultra Phenyl 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9105518	Ultra Phenyl 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9105519	Ultra Phenyl 5um	100 x 30mm Confirm customer has evaluated an analytical size column
9105521	Ultra Phenyl 5um	200 x 1.0mm
9105522	Ultra Phenyl 5um	200 x 2.1mm
9105522-700	Ultra Phenyl 5um	200 x 2.1mm with Trident
9105523	Ultra Phenyl 5um	200 x 3.2mm
9105523-700	Ultra Phenyl 5um	200 x 3.2mm with Trident
9105525	Ultra Phenyl 5um	200 x 4.6mm
9105525-700	Ultra Phenyl 5um	200 x 4.6mm with Trident
9105531	Ultra Phenyl 5um	30 x 1.0mm
9105532	Ultra Phenyl 5um	30 x 2.1mm
9105532-700	Ultra Phenyl 5um	30 x 2.1mm with Trident
9105533	Ultra Phenyl 5um	30 x 3.2mm
9105533-700	Ultra Phenyl 5um	30 x 3.2mm with Trident
9105535	Ultra Phenyl 5um	30 x 4.6mm
9105535-700	Ultra Phenyl 5um	30 x 4.6mm with Trident
9105550	Ultra Phenyl 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9105551	Ultra Phenyl 5um	50 x 1.0mm
9105552	Ultra Phenyl 5um	50 x 2.1mm
9105552-700	Ultra Phenyl 5um	50 x 2.1mm with Trident
9105553	Ultra Phenyl 5um	50 x 3.2mm
9105553-700	Ultra Phenyl 5um	50 x 3.2mm with Trident
9105555	Ultra Phenyl 5um	50 x 4.6mm
9105555-700	Ultra Phenyl 5um	50 x 4.6mm with Trident
9105557	Ultra Phenyl 5um	50 x 10mm Confirm customer has evaluated an analytical size column
9105558	Ultra Phenyl 5um	50 x 21.2mm Confirm customer has evaluated an analytical size column
9105559	Ultra Phenyl 5um	50 x 30mm Confirm customer has evaluated an analytical size column
9105560	Ultra Phenyl 5um	150 x 50mm Confirm customer has evaluated an analytical size column
9105561	Ultra Phenyl 5um	150 x 1.0mm
9105562	Ultra Phenyl 5um	150 x 2.1mm
9105562-700	Ultra Phenyl 5um	150 x 2.1mm with Trident
9105563	Ultra Phenyl 5um	150 x 3.2mm
9105563-700	Ultra Phenyl 5um	150 x 3.2mm with Trident
9105565	Ultra Phenyl 5um	150 x 4.6mm
9105565-700	Ultra Phenyl 5um	150 x 4.6mm with Trident
9105567	Ultra Phenyl 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9105568	Ultra Phenyl 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9105569	Ultra Phenyl 5um	150 x 30mm Confirm customer has evaluated an analytical size column
9105570	Ultra Phenyl 5um	250 x 50mm Confirm customer has evaluated an analytical size column
9105571	Ultra Phenyl 5um	250 x 1.0mm
9105572	Ultra Phenyl 5um	250 x 2.1mm
9105572-700	Ultra Phenyl 5um	250 x 2.1mm with Trident
9105573	Ultra Phenyl 5um	250 x 3.2mm
9105573-700	Ultra Phenyl 5um	250 x 3.2mm with Trident
9105575	Ultra Phenyl 5um	250 x 4.6mm
9105575-700	Ultra Phenyl 5um	250 x 4.6mm with Trident
9105577	Ultra Phenyl 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9105578	Ultra Phenyl 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column
9105579	Ultra Phenyl 5um	250 x 30mm Confirm customer has evaluated an analytical size column
9106311	Ultra Cyano 3um	100 x 1.0mm
9106312	Ultra Cyano 3um	100 x 2.1mm
9106312-700	Ultra Cyano 3um	100 x 2.1mm with Trident
9106313	Ultra Cyano 3um	100 x 3.2mm

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

 HROMalytic +61(0)3 9762 2034 ECHnology Pty Ltd		Australian Distributors Importers & Manufacturers www.chromtech.net.au
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA		

(Ctrl F)
> Cat.#s

> Current System AUD Prices

