

9160579	Allure Silica 5um	250 x 30mm Confirm customer has evaluated an analytical size column
91615	Allure Basix 5um	Bulk Packing Material
91615	Allure Basix 5um	Bulk Packing Material
91615	Allure Basix 5um	Bulk Packing Material
91615	Allure Basix 5um	Bulk Packing Material
916150210	Allure Basix	10 x 4.0mm 3-pk
916150212	Allure Basix	10 x 2.1mm 3-pk
916150220	Allure Basix	20 x 4.0mm 2-pk
916150222	Allure Basix	20 x 2.1mm 2-pk
9161510	Allure Basix 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9161511	Allure Basix 5um	100 x 1.0mm
9161512	Allure Basix 5um	100 x 2.1mm
9161512-700	Allure Basix 5um	100 x 2.1mm with Trident
9161513	Allure Basix 5um	100 x 3.2mm
9161513-700	Allure Basix 5um	100 x 3.2mm with Trident
9161515	Allure Basix 5um	100 x 4.6mm
9161515-700	Allure Basix 5um	100 x 4.6mm with Trident
9161517	Allure Basix 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9161518	Allure Basix 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9161519	Allure Basix 5um	100 x 30mm Confirm customer has evaluated an analytical size column
9161521	Allure Basix 5um	200 x 1.0mm
9161522	Allure Basix 5um	200 x 2.1mm
9161522-700	Allure Basix 5um	200 x 2.1mm with Trident
9161523	Allure Basix 5um	200 x 3.2mm
9161523-700	Allure Basix 5um	200 x 3.2mm with Trident
9161525	Allure Basix 5um	200 x 4.6mm
9161525-700	Allure Basix 5um	200 x 4.6mm with Trident
9161531	Allure Basix 5um	30 x 1.0mm
9161532	Allure Basix 5um	30 x 2.1mm
9161532-700	Allure Basix 5um	30 x 2.1mm with Trident
9161533	Allure Basix 5um	30 x 3.2mm
9161533-700	Allure Basix 5um	30 x 3.2mm with Trident
9161535	Allure Basix 5um	30 x 4.6mm
9161535-700	Allure Basix 5um	30 x 4.6mm with Trident
9161550	Allure Basix 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9161551	Allure Basix 5um	50 x 1.0mm
9161552	Allure Basix 5um	50 x 2.1mm
9161552-700	Allure Basix 5um	50 x 2.1mm with Trident
9161553	Allure Basix 5um	50 x 3.2mm
9161553-700	Allure Basix 5um	50 x 3.2mm with Trident
9161555	Allure Basix 5um	50 x 4.6mm
9161555-700	Allure Basix 5um	50 x 4.6mm with Trident
9161557	Allure Basix 5um	50 x 10mm Confirm customer has evaluated an analytical size column
9161558	Allure Basix 5um	50 x 21.2mm Confirm customer has evaluated an analytical size column
9161559	Allure Basix 5um	50 x 30mm Confirm customer has evaluated an analytical size column
9161560	Allure Basix 5um	150 x 50mm Confirm customer has evaluated an analytical size column
9161561	Allure Basix 5um	150 x 1.0mm
9161562	Allure Basix 5um	150 x 2.1mm
9161562-700	Allure Basix 5um	150 x 2.1mm with Trident
9161563	Allure Basix 5um	150 x 3.2mm
9161563-700	Allure Basix 5um	150 x 3.2mm with Trident
9161565	Allure Basix 5um	150 x 4.6mm
9161565-700	Allure Basix 5um	150 x 4.6mm with Trident
9161567	Allure Basix 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9161568	Allure Basix 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9161569	Allure Basix 5um	150 x 30mm Confirm customer has evaluated an analytical size column
9161570	Allure Basix 5um	250 x 50mm Confirm customer has evaluated an analytical size column
9161571	Allure Basix 5um	250 x 1.0mm
9161572	Allure Basix 5um	250 x 2.1mm
9161572-700	Allure Basix 5um	250 x 2.1mm with Trident
9161573	Allure Basix 5um	250 x 3.2mm
9161573-700	Allure Basix 5um	250 x 3.2mm with Trident
9161575	Allure Basix 5um	250 x 4.6mm
9161575-700	Allure Basix 5um	250 x 4.6mm with Trident
9161577	Allure Basix 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9161578	Allure Basix 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
 Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

(Ctrl F)
 > Cat.#s

> Current System AUD Prices

