

9166568	Allure BiPh 5um	150 x 21.2mm Confirm customer has
9166569	Allure BiPh 5um	150 x 30mm Confirm customer has ev
9166570	Allure BiPh 5um	250 x 50mm Confirm customer has ev
9166571	Allure BiPh 5um	250 x 1.0mm
9166572	Allure BiPh 5um	250 x 2.1mm
9166572-700	Allure BiPh 5um	250 x 2.1mm with Trident
9166573	Allure BiPh 5um	250 x 3.2mm
9166573-700	Allure BiPh 5um	250 x 3.2mm with Trident
9166575	Allure BiPh 5um	250 x 4.6mm
9166575-700	Allure BiPh 5um	250 x 4.6mm with Trident
9166577	Allure BiPh 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9166578	Allure BiPh 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column
9166579	Allure BiPh 5um	250 x 30mm Confirm customer has evaluated an analytical size column
916850210	Allure AQ C18	10 x 4.0mm 3-pk
916850212	Allure AQ C18	10 x 2.1mm 3-pk
916850220	Allure AQ C18	20 x 4.0mm 2-pk
916850222	Allure AQ C18	20 x 2.1mm 2-pk
9168511	Allure AQ C18 5um	100 x 1.0mm
9168512	Allure AQ C18 5um	100 x 2.1mm
9168512-700	Allure AQ C18 5um	100 x 2.1mm with Trident
9168513	Allure AQ C18 5um	100 x 3.2mm
9168513-700	Allure AQ C18 5um	100 x 3.2mm with Trident
9168515	Allure AQ C18 5um	100 x 4.6mm
9168515-700	Allure AQ C18 5um	100 x 4.6mm with Trident
9168521	Allure AQ C18 5um	200 x 1.0mm
9168522	Allure AQ C18 5um	200 x 2.1mm
9168522-700	Allure AQ C18 5um	200 x 2.1mm with Trident
9168523	Allure AQ C18 5um	200 x 3.2mm
9168523-700	Allure AQ C18 5um	200 x 3.2mm with Trident
9168525	Allure AQ C18 5um	200 x 4.6mm
9168525-700	Allure AQ C18 5um	200 x 4.6mm with Trident
9168531	Allure AQ C18 5um	30 x 1.0mm
9168532	Allure AQ C18 5um	30 x 2.1mm
9168532-700	Allure AQ C18 5um	30 x 2.1mm with Trident
9168533	Allure AQ C18 5um	30 x 3.2mm
9168533-700	Allure AQ C18 5um	30 x 3.2mm with Trident
9168535	Allure AQ C18 5um	30 x 4.6mm
9168535-700	Allure AQ C18 5um	30 x 4.6mm with Trident
9168551	Allure AQ C18 5um	50 x 1.0mm
9168552	Allure AQ C18 5um	50 x 2.1mm
9168552-700	Allure AQ C18 5um	50 x 2.1mm with Trident
9168553	Allure AQ C18 5um	50 x 3.2mm
9168553-700	Allure AQ C18 5um	50 x 3.2mm with Trident
9168555	Allure AQ C18 5um	50 x 4.6mm
9168555-700	Allure AQ C18 5um	50 x 4.6mm with Trident
9168561	Allure AQ C18 5um	150 x 1.0mm
9168562	Allure AQ C18 5um	150 x 2.1mm
9168562-700	Allure AQ C18 5um	150 x 2.1mm with Trident
9168563	Allure AQ C18 5um	150 x 3.2mm
9168563-700	Allure AQ C18 5um	150 x 3.2mm with Trident
9168565	Allure AQ C18 5um	150 x 4.6mm
9168565-700	Allure AQ C18 5um	150 x 4.6mm with Trident
9168571	Allure AQ C18 5um	250 x 1.0mm
9168572	Allure AQ C18 5um	250 x 2.1mm
9168572-700	Allure AQ C18 5um	250 x 2.1mm with Trident
9168573	Allure AQ C18 5um	250 x 3.2mm
9168573-700	Allure AQ C18 5um	250 x 3.2mm with Trident
9168575	Allure AQ C18 5um	250 x 4.6mm
9168575-700	Allure AQ C18 5um	250 x 4.6mm with Trident
916950210	Allure PFPP	10 x 4.0mm 3-pk
916950212	Allure PFPP	10 x 2.1mm 3-pk
916950220	Allure PFPP	20 x 4.0mm 2-pk
916950222	Allure PFPP	20 x 2.1mm 2-pk
9169510	Allure PFPP 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9169511	Allure PFPP 5um	100 x 1.0mm
9169512	Allure PFPP 5um	100 x 2.1mm

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

 HROMalytic +61(0)3 9762 2034 ECHnology Pty Ltd		Australian Distributors Importers & Manufacturers www.chromtech.net.au
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA		

(Ctrl F)
> Cat.#s

> Current System AUD Prices

