

9174332	Ultra C18 3um	30 x 2.1mm
9174332-700	Ultra C18 3um	30 x 2.1mm with Trident
9174333	Ultra C18 3um	30 x 3.2mm
9174333-700	Ultra C18 3um	30 x 3.2mm with Trident
9174335	Ultra C18 3um	30 x 4.6mm
9174335-700	Ultra C18 3um	30 x 4.6mm with Trident
9174351	Ultra C18 3um	50 x 1.0mm
9174352	Ultra C18 3um	50 x 2.1mm
9174352-700	Ultra C18 3um	50 x 2.1mm with Trident
9174353	Ultra C18 3um	50 x 3.2mm
9174353-700	Ultra C18 3um	50 x 3.2mm with Trident
9174355	Ultra C18 3um	50 x 4.6mm
9174355-700	Ultra C18 3um	50 x 4.6mm with Trident
917435B	Ultra C18 3um	50 x 0.3mm Capillary Column
9174361	Ultra C18 3um	150 x 1.0mm
9174362	Ultra C18 3um	150 x 2.1mm
9174362-700	Ultra C18 3um	150 x 2.1mm with Trident
9174363	Ultra C18 3um	150 x 3.2mm
9174363-700	Ultra C18 3um	150 x 3.2mm with Trident
9174365	Ultra C18 3um	150 x 4.6mm
9174365-700	Ultra C18 3um	150 x 4.6mm with Trident
917436B	Ultra C18 3um	150 x 0.3mm Capillary Column
91745	Ultra C18 5um	Bulk Packing Material 10 grams
91745	Ultra C18 5um	Bulk Packing Material 10 grams
91745	Ultra C18 5um	Bulk Packing Material 10 grams
917450210	Ultra C18	10 x 4.0mm 3-pk
917450212	Ultra C18	10 x 2.1mm 3-pk
917450220	Ultra C18	20 x 4.0mm 2-pk
917450222	Ultra C18	20 x 2.1mm 2-pk
9174510	Ultra C18 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9174511	Ultra C18 5um	100 x 1.0mm
9174512	Ultra C18 5um	100 x 2.1mm
9174512-700	Ultra C18 5um	100 x 2.1mm with Trident
9174513	Ultra C18 5um	100 x 3.2mm
9174513-700	Ultra C18 5um	100 x 3.2mm with Trident
9174514	Ultra C18 5um	100 x 4.0mm
9174514-700	Ultra C18 5um	100 x 4.0mm with Trident
9174515	Ultra C18 5um	100 x 4.6mm
9174515-700	Ultra C18 5um	100 x 4.6mm with Trident
9174517	Ultra C18 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9174518	Ultra C18 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9174519	Ultra C18 5um	100 x 30mm Confirm customer has evaluated an analytical size column
9174521	Ultra C18 5um	200 x 1.0mm
9174522	Ultra C18 5um	200 x 2.1mm
9174522-700	Ultra C18 5um	200 x 2.1mm with Trident
9174523	Ultra C18 5um	200 x 3.2mm
9174523-700	Ultra C18 5um	200 x 3.2mm with Trident
9174525	Ultra C18 5um	200 x 4.6mm
9174525-700	Ultra C18 5um	200 x 4.6mm with Trident
9174531	Ultra C18 5um	30 x 1.0mm
9174532	Ultra C18 5um	30 x 2.1mm
9174532-700	Ultra C18 5um	30 x 2.1mm with Trident
9174533	Ultra C18 5um	30 x 3.2mm
9174533-700	Ultra C18 5um	30 x 3.2mm with Trident
9174535	Ultra C18 5um	30 x 4.6mm
9174535-700	Ultra C18 5um	30 x 4.6mm with Trident
9174550	Ultra C18 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9174551	Ultra C18 5um	50 x 1.0mm
9174552	Ultra C18 5um	50 x 2.1mm
9174552-700	Ultra C18 5um	50 x 2.1mm with Trident
9174553	Ultra C18 5um	50 x 3.2mm
9174553-700	Ultra C18 5um	50 x 3.2mm with Trident
9174555	Ultra C18 5um	50 x 4.6mm
9174555-700	Ultra C18 5um	50 x 4.6mm with Trident
9174557	Ultra C18 5um	50 x 10mm Confirm customer has evaluated an analytical size column

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
 Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

(Ctrl F)
> Cat.#s

> Current System AUD Prices

