

9176563-700	Ultra PFP 5um	150 x 3.2mm with Trident
9176565	Ultra PFP 5um	150 x 4.6mm
9176565-700	Ultra PFP 5um	150 x 4.6mm with Trident
9176567	Ultra PFP 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9176568	Ultra PFP 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9176569	Ultra PFP 5um	150 x 30mm Confirm customer has evaluated an analytical size column
9176570	Ultra PFP 5um	250 x 50mm Confirm customer has evaluated an analytical size column
9176571	Ultra PFP 5um	250 x 1.0mm
9176572	Ultra PFP 5um	250 x 2.1mm
9176572-700	Ultra PFP 5um	250 x 2.1mm with Trident
9176573	Ultra PFP 5um	250 x 3.2mm
9176573-700	Ultra PFP 5um	250 x 3.2mm with Trident
9176575	Ultra PFP 5um	250 x 4.6mm
9176575-700	Ultra PFP 5um	250 x 4.6mm with Trident
9176577	Ultra PFP 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9176578	Ultra PFP 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column
9176579	Ultra PFP 5um	250 x 30mm Confirm customer has evaluated an analytical size column
9177311	Ultra Carbamate 3um	100 x 1.0mm
9177312	Ultra Carbamate 3um	100 x 2.1mm
9177313	Ultra Carbamate 3um	100 x 3.2mm
9177315	Ultra Carbamate 3um	100 x 4.6mm
9177315-700	Ultra Carbamate 3um	100 x 4.6mm with Trident
9177351	Ultra Carbamate 3um	50 x 1.0mm
9177352	Ultra Carbamate 3um	50 x 2.1mm
9177353	Ultra Carbamate 3um	50 x 3.2mm
9177354	Ultra Carbamate 3um	50 x 4.0mm
9177355	Ultra Carbamate 3um	50 x 4.6mm
9177355-700	Ultra Carbamate 3um	50 x 4.6mm with Trident
917750210	Ultra Carbamate	10 x 4.0mm 3-pk
917750212	Ultra Carbamate	10 x 2.1mm 3-pk
917750220	Ultra Carbamate	20 x 4.0mm 2-pk
917750222	Ultra Carbamate	20 x 2.1mm 2-pk
9177575	Ultra Carbamate 5um	250 x 4.6mm
9177575-700	Ultra Carbamate 5um	250 x 4.6mm with Trident
9178311	Ultra AQ C18 3um	100 x 1.0mm
9178312	Ultra AQ C18 3um	100 x 2.1mm
9178312-700	Ultra AQ C18 3um	100 x 2.1mm with Trident
9178313	Ultra AQ C18 3um	100 x 3.2mm
9178313-700	Ultra AQ C18 3um	100 x 3.2mm with Trident
9178315	Ultra AQ C18 3um	100 x 4.6mm
9178315-700	Ultra AQ C18 3um	100 x 4.6mm with Trident
9178331	Ultra AQ C18 3um	30 x 1.0mm
9178332	Ultra AQ C18 3um	30 x 2.1mm
9178332-700	Ultra AQ C18 3um	30 x 2.1mm with Trident
9178333	Ultra AQ C18 3um	30 x 3.2mm
9178333-700	Ultra AQ C18 3um	30 x 3.2mm with Trident
9178335	Ultra AQ C18 3um	30 x 4.6mm
9178335-700	Ultra AQ C18 3um	30 x 4.6mm with Trident
9178351	Ultra AQ C18 3um	50 x 1.0mm
9178352	Ultra AQ C18 3um	50 x 2.1mm
9178352-700	Ultra AQ C18 3um	50 x 2.1mm with Trident
9178353	Ultra AQ C18 3um	50 x 3.2mm
9178353-700	Ultra AQ C18 3um	50 x 3.2mm with Trident
9178355	Ultra AQ C18 3um	50 x 4.6mm
9178355-700	Ultra AQ C18 3um	50 x 4.6mm with Trident
917835B	Ultra AQ C18 3um	50 x 0.3mm Capillary Column
9178361	Ultra AQ C18 3um	150 x 1.0mm
9178362	Ultra AQ C18 3um	150 x 2.1mm
9178362-700	Ultra AQ C18 3um	150 x 2.1mm with Trident
9178363	Ultra AQ C18 3um	150 x 3.2mm
9178363-700	Ultra AQ C18 3um	150 x 3.2mm with Trident
9178365	Ultra AQ C18 3um	150 x 4.6mm
9178365-700	Ultra AQ C18 3um	150 x 4.6mm with Trident
917836B	Ultra AQ C18 3um	150 x 0.3mm Capillary Column
917850210	Ultra AQ C18	10 x 4.0mm 3-pk
917850212	Ultra AQ C18	10 x 2.1mm 3-pk

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

		<b>Australian Distributors</b> Importers & Manufacturers <a href="http://www.chromtech.net.au">www.chromtech.net.au</a>
Website NEW : <a href="http://www.chromalytic.net.au">www.chromalytic.net.au</a> E-mail : <a href="mailto:info@chromtech.net.au">info@chromtech.net.au</a> Tel: 03 9762 2034 . . . in AUSTRALIA		

( Ctrl F )  
> Cat.#s

> Current System AUD Prices

