

917850220	Ultra AQ C18	20 x 4.0mm 2-pk
917850222	Ultra AQ C18	20 x 2.1mm 2-pk
9178510	Ultra AQ C18 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9178511	Ultra AQ C18 5um	100 x 1.0mm
9178512	Ultra AQ C18 5um	100 x 2.1mm
9178512-700	Ultra AQ C18 5um	100 x 2.1mm with Trident
9178513	Ultra AQ C18 5um	100 x 3.2mm
9178513-700	Ultra AQ C18 5um	100 x 3.2mm with Trident
9178515	Ultra AQ C18 5um	100 x 4.6mm
9178515-700	Ultra AQ C18 5um	100 x 4.6mm with Trident
9178517	Ultra AQ C18 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9178518	Ultra AQ C18 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9178519	Ultra AQ C18 5um	100 x 30mm Confirm customer has evaluated an analytical size column
9178521	Ultra AQ C18 5um	200 x 1.0mm
9178522	Ultra AQ C18 5um	200 x 2.1mm
9178522-700	Ultra AQ C18 5um	200 x 2.1mm with Trident
9178523	Ultra AQ C18 5um	200 x 3.2mm
9178523-700	Ultra AQ C18 5um	200 x 3.2mm with Trident
9178525	Ultra AQ C18 5um	200 x 4.6mm
9178525-700	Ultra AQ C18 5um	200 x 4.6mm with Trident
9178531	Ultra AQ C18 5um	30 x 1.0mm
9178532	Ultra AQ C18 5um	30 x 2.1mm
9178532-700	Ultra AQ C18 5um	30 x 2.1mm with Trident
9178533	Ultra AQ C18 5um	30 x 3.2mm
9178533-700	Ultra AQ C18 5um	30 x 3.2mm with Trident
9178535	Ultra AQ C18 5um	30 x 4.6mm
9178535-700	Ultra AQ C18 5um	30 x 4.6mm with Trident
9178550	Ultra AQ C18 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9178551	Ultra AQ C18 5um	50 x 1.0mm
9178552	Ultra AQ C18 5um	50 x 2.1mm
9178552-700	Ultra AQ C18 5um	50 x 2.1mm with Trident
9178553	Ultra AQ C18 5um	50 x 3.2mm
9178553-700	Ultra AQ C18 5um	50 x 3.2mm with Trident
9178555	Ultra AQ C18 5um	50 x 4.6mm
9178555-700	Ultra AQ C18 5um	50 x 4.6mm with Trident
9178557	Ultra AQ C18 5um	50 x 10mm Confirm customer has evaluated an analytical size column
9178558	Ultra AQ C18 5um	50 x 21.2mm Confirm customer has evaluated an analytical size column
9178559	Ultra AQ C18 5um	50 x 30mm Confirm customer has evaluated an analytical size column
9178560	Ultra AQ C18 5um	150 x 50mm Confirm customer has evaluated an analytical size column
9178561	Ultra AQ C18 5um	150 x 1.0mm
9178562	Ultra AQ C18 5um	150 x 2.1mm
9178562-700	Ultra AQ C18 5um	150 x 2.1mm with Trident
9178563	Ultra AQ C18 5um	150 x 3.2mm
9178563-700	Ultra AQ C18 5um	150 x 3.2mm with Trident
9178565	Ultra AQ C18 5um	150 x 4.6mm
9178565-700	Ultra AQ C18 5um	150 x 4.6mm with Trident
9178567	Ultra AQ C18 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9178568	Ultra AQ C18 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9178569	Ultra AQ C18 5um	150 x 30mm Confirm customer has evaluated an analytical size column
9178570	Ultra AQ C18 5um	250 x 50mm Confirm customer has evaluated an analytical size column
9178571	Ultra AQ C18 5um	250 x 1.0mm
9178572	Ultra AQ C18 5um	250 x 2.1mm
9178572-700	Ultra AQ C18 5um	250 x 2.1mm with Trident
9178573	Ultra AQ C18 5um	250 x 3.2mm
9178573-700	Ultra AQ C18 5um	250 x 3.2mm with Trident
9178575	Ultra AQ C18 5um	250 x 4.6mm
9178575-700	Ultra AQ C18 5um	250 x 4.6mm with Trident
9178577	Ultra AQ C18 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9178578	Ultra AQ C18 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column
9178579	Ultra AQ C18 5um	250 x 30mm Confirm customer has evaluated an analytical size column
9178584	Ultra AQ C18 5um	300 x 4.0mm
9179311	Ultra PFPP 3um	100 x 1.0mm
9179312	Ultra PFPP 3um	100 x 2.1mm
9179312-700	Ultra PFPP 3um	100 x 2.1mm with Trident
9179313	Ultra PFPP 3um	100 x 3.2mm
9179313-700	Ultra PFPP 3um	100 x 3.2mm with Trident

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

		Australian Distributors Importers & Manufacturers www.chromtech.net.au
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA		

(Ctrl F)
> Cat.#s

> Current System AUD Prices

