

92105	Pinnacle II Silica 5um	Bulk Packing Material
92105	Pinnacle II Silica 5um	Bulk Packing Material
921050210	Pinnacle II Silica	10 x 4.0mm 3-pk
921050212	Pinnacle II Silica	10 x 2.1mm 3-pk
921050220	Pinnacle II Silica	20 x 4.0mm 2-pk
921050222	Pinnacle II Silica	20 x 2.1mm 2-pk
9210510	Pinnacle II Silica 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9210511	Pinnacle II Silica 5um	100 x 1.0mm
9210512	Pinnacle II Silica 5um	100 x 2.1mm
9210512-700	Pinnacle II Silica 5um	100 x 2.1mm with Trident
9210513	Pinnacle II Silica 5um	100 x 3.2mm
9210513-700	Pinnacle II Silica 5um	100 x 3.2mm with Trident
9210515	Pinnacle II Silica 5um	100 x 4.6mm
9210515-700	Pinnacle II Silica 5um	100 x 4.6mm with Trident
9210517	Pinnacle II Silica 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9210518	Pinnacle II Silica 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9210519	Pinnacle II Silica 5um	100 x 30mm Confirm customer has evaluated an analytical size column
9210521	Pinnacle II Silica 5um	200 x 1.0mm
9210522	Pinnacle II Silica 5um	200 x 2.1mm
9210522-700	Pinnacle II Silica 5um	200 x 2.1mm with Trident
9210523	Pinnacle II Silica 5um	200 x 3.2mm
9210523-700	Pinnacle II Silica 5um	200 x 3.2mm with Trident
9210525	Pinnacle II Silica 5um	200 x 4.6mm
9210525-700	Pinnacle II Silica 5um	200 x 4.6mm with Trident
9210531	Pinnacle II Silica 5um	30 x 1.0mm
9210532	Pinnacle II Silica 5um	30 x 2.1mm
9210532-700	Pinnacle II Silica 5um	30 x 2.1mm with Trident
9210533	Pinnacle II Silica 5um	30 x 3.2mm
9210533-700	Pinnacle II Silica 5um	30 x 3.2mm with Trident
9210535	Pinnacle II Silica 5um	30 x 4.6mm
9210535-700	Pinnacle II Silica 5um	30 x 4.6mm with Trident
9210550	Pinnacle II Silica 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9210551	Pinnacle II Silica 5um	50 x 1.0mm
9210552	Pinnacle II Silica 5um	50 x 2.1mm
9210552-700	Pinnacle II Silica 5um	50 x 2.1mm with Trident
9210553	Pinnacle II Silica 5um	50 x 3.2mm
9210553-700	Pinnacle II Silica 5um	50 x 3.2mm with Trident
9210555	Pinnacle II Silica 5um	50 x 4.6mm
9210555-700	Pinnacle II Silica 5um	50 x 4.6mm with Trident
9210557	Pinnacle II Silica 5um	50 x 10mm Confirm customer has evaluated an analytical size column
9210558	Pinnacle II Silica 5um	50 x 21.2mm Confirm customer has evaluated an analytical size column
9210559	Pinnacle II Silica 5um	50 x 30mm Confirm customer has evaluated an analytical size column
9210560	Pinnacle II Silica 5um	150 x 50mm Confirm customer has evaluated an analytical size column
9210561	Pinnacle II Silica 5um	150 x 1.0mm
9210562	Pinnacle II Silica 5um	150 x 2.1mm
9210562-700	Pinnacle II Silica 5um	150 x 2.1mm with Trident
9210563	Pinnacle II Silica 5um	150 x 3.2mm
9210563-700	Pinnacle II Silica 5um	150 x 3.2mm with Trident
9210565	Pinnacle II Silica 5um	150 x 4.6mm
9210565-700	Pinnacle II Silica 5um	150 x 4.6mm with Trident
9210567	Pinnacle II Silica 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9210568	Pinnacle II Silica 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9210569	Pinnacle II Silica 5um	150 x 30mm Confirm customer has evaluated an analytical size column
9210570	Pinnacle II Silica 5um	250 x 50mm Confirm customer has evaluated an analytical size column
9210571	Pinnacle II Silica 5um	250 x 1.0mm
9210572	Pinnacle II Silica 5um	250 x 2.1mm
9210572-700	Pinnacle II Silica 5um	250 x 2.1mm with Trident
9210573	Pinnacle II Silica 5um	250 x 3.2mm
9210573-700	Pinnacle II Silica 5um	250 x 3.2mm with Trident
9210575	Pinnacle II Silica 5um	250 x 4.6mm
9210575-700	Pinnacle II Silica 5um	250 x 4.6mm with Trident
9210577	Pinnacle II Silica 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9210578	Pinnacle II Silica 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column
9210579	Pinnacle II Silica 5um	250 x 30mm Confirm customer has evaluated an analytical size column
9212078	Pinnacle II C4 10um	250 x 21.2mm Confirm customer has evaluated an analytical size column
92133	Pinnacle II C8 3um	Bulk Packing Material

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
 Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

(Ctrl F)
> Cat.#s

> Current System AUD Prices

