

9214363-700	Pinnacle II C18 3um	150 x 3.2mm with Trident
9214365	Pinnacle II C18 3um	150 x 4.6mm
9214365-700	Pinnacle II C18 3um	150 x 4.6mm with Trident
92145	Pinnacle II C18 5um	Bulk Packing Material
92145	Pinnacle II C18 5um	Bulk Packing Material
92145	Pinnacle II C18 5um	Bulk Packing Material
921450210	Pinnacle II C18	10 x 4.0mm 3-pk
921450212	Pinnacle II C18	10 x 2.1mm 3-pk
921450220	Pinnacle II C18	20 x 4.0mm 2-pk
921450222	Pinnacle II C18	20 x 2.1mm 2-pk
9214510	Pinnacle II C18 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9214511	Pinnacle II C18 5um	100 x 1.0mm
9214512	Pinnacle II C18 5um	100 x 2.1mm
9214512-700	Pinnacle II C18 5um	100 x 2.1mm with Trident
9214513	Pinnacle II C18 5um	100 x 3.2mm
9214513-700	Pinnacle II C18 5um	100 x 3.2mm with Trident
9214514	Pinnacle II C18 5um	100 x 4.0mm
9214515	Pinnacle II C18 5um	100 x 4.6mm
9214515-700	Pinnacle II C18 5um	100 x 4.6mm with Trident
9214517	Pinnacle II C18 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9214518	Pinnacle II C18 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9214519	Pinnacle II C18 5um	100 x 30mm Confirm customer has evaluated an analytical size column
9214521	Pinnacle II C18 5um	200 x 1.0mm
9214522	Pinnacle II C18 5um	200 x 2.1mm
9214522-700	Pinnacle II C18 5um	200 x 2.1mm with Trident
9214523	Pinnacle II C18 5um	200 x 3.2mm
9214523-700	Pinnacle II C18 5um	200 x 3.2mm with Trident
9214525	Pinnacle II C18 5um	200 x 4.6mm
9214525-700	Pinnacle II C18 5um	200 x 4.6mm with Trident
9214531	Pinnacle II C18 5um	30 x 1.0mm
9214532	Pinnacle II C18 5um	30 x 2.1mm
9214532-700	Pinnacle II C18 5um	30 x 2.1mm with Trident
9214533	Pinnacle II C18 5um	30 x 3.2mm
9214533-700	Pinnacle II C18 5um	30 x 3.2mm with Trident
9214535	Pinnacle II C18 5um	30 x 4.6mm
9214535-700	Pinnacle II C18 5um	30 x 4.6mm with Trident
9214550	Pinnacle II C18 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9214551	Pinnacle II C18 5um	50 x 1.0mm
9214552	Pinnacle II C18 5um	50 x 2.1mm
9214552-700	Pinnacle II C18 5um	50 x 2.1mm with Trident
9214553	Pinnacle II C18 5um	50 x 3.2mm
9214553-700	Pinnacle II C18 5um	50 x 3.2mm with Trident
9214555	Pinnacle II C18 5um	50 x 4.6mm
9214555-700	Pinnacle II C18 5um	50 x 4.6mm with Trident
9214557	Pinnacle II C18 5um	50 x 10mm Confirm customer has evaluated an analytical size column
9214558	Pinnacle II C18 5um	50 x 21.2mm Confirm customer has evaluated an analytical size column
9214559	Pinnacle II C18 5um	50 x 30mm Confirm customer has evaluated an analytical size column
9214560	Pinnacle II C18 5um	150 x 50mm Confirm customer has evaluated an analytical size column
9214561	Pinnacle II C18 5um	150 x 1.0mm
9214562	Pinnacle II C18 5um	150 x 2.1mm
9214562-700	Pinnacle II C18 5um	150 x 2.1mm with Trident
9214563	Pinnacle II C18 5um	150 x 3.2mm
9214563-700	Pinnacle II C18 5um	150 x 3.2mm with Trident
9214564	Pinnacle II C18 5um	150 x 4.0mm
9214564-700	Pinnacle II C18 5um	150 x 4.0mm with Trident
9214565	Pinnacle II C18 5um	150 x 4.6mm
9214565-700	Pinnacle II C18 5um	150 x 4.6mm with Trident
9214567	Pinnacle II C18 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9214568	Pinnacle II C18 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9214569	Pinnacle II C18 5um	150 x 30mm Confirm customer has evaluated an analytical size column
9214570	Pinnacle II C18 5um	250 x 50mm Confirm customer has evaluated an analytical size column
9214571	Pinnacle II C18 5um	250 x 1.0mm
9214572	Pinnacle II C18 5um	250 x 2.1mm
9214572-700	Pinnacle II C18 5um	250 x 2.1mm with Trident
9214573	Pinnacle II C18 5um	250 x 3.2mm

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
 Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

(Ctrl F)
> Cat.#s

> Current System AUD Prices

