

9603570	Ultra II C8 5um	250 x 50mm Confirm customer has ev
9603571	Ultra II C8 5um	250 x 1.0mm
9603572	Ultra II C8 5um	250 x 2.1mm
9603572-700	Ultra II C8 5um	250 x 2.1mm With Trident
9603573	Ultra II C8 5um	250 x 3.2mm
9603573-700	Ultra II C8 5um	250 x 3.2mm With Trident
9603575	Ultra II C8 5um	250 x 4.6mm
9603575-700	Ultra II C8 5um	250 x 4.6mm With Trident
9603577	Ultra II C8 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9603578	Ultra II C8 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column
9603579	Ultra II C8 5um	250 x 30mm Confirm customer has evaluated an analytical size column
960357E	Ultra II C8 5um	250 x 3.0mm
9604077	Ultra II C18 10um	250 x 10mm Confirm customer has evaluated an analytical size column
9604212	Ultra II C18 1.9um	100 x 2.1mm
960421E	Ultra II C18 1.9um	100 x 3.0mm
9604232	Ultra II C18 1.9um	30 x 2.1mm
960423E	Ultra II C18 1.9um	30 x 3.0mm
9604252	Ultra II C18 1.9um	50 x 2.1mm
960425E	Ultra II C18 1.9um	50 x 3.0mm
96043	Ultra II C18 3um	Bulk Packing Material
96043	Ultra II C18 3um	Bulk Packing Material
96043	Ultra II C18 3um	Bulk Packing Material
9604311	Ultra II C18 3um	100 x 1.0mm
9604312	Ultra II C18 3um	100 x 2.1mm
9604312-700	Ultra II C18 3um	100 x 2.1mm With Trident
9604313	Ultra II C18 3um	100 x 3.2mm
9604313-700	Ultra II C18 3um	100 x 3.2mm With Trident
9604315	Ultra II C18 3um	100 x 4.6mm
9604315-700	Ultra II C18 3um	100 x 4.6mm With Trident
960431E	Ultra II C18 3um	100 x 3.0mm
9604331	Ultra II C18 3um	30 x 1.0mm
9604332	Ultra II C18 3um	30 x 2.1mm
9604332-700	Ultra II C18 3um	30 x 2.1mm With Trident
9604333	Ultra II C18 3um	30 x 3.2mm
9604333-700	Ultra II C18 3um	30 x 3.2mm With Trident
9604335	Ultra II C18 3um	30 x 4.6mm
9604335-700	Ultra II C18 3um	30 x 4.6mm With Trident
960433E	Ultra II C18 3um	30 x 3.0mm
9604351	Ultra II C18 3um	50 x 1.0mm
9604352	Ultra II C18 3um	50 x 2.1mm
9604352-700	Ultra II C18 3um	50 x 2.1mm With Trident
9604353	Ultra II C18 3um	50 x 3.2mm
9604353-700	Ultra II C18 3um	50 x 3.2mm With Trident
9604355	Ultra II C18 3um	50 x 4.6mm
9604355-700	Ultra II C18 3um	50 x 4.6mm With Trident
960435E	Ultra II C18 3um	50 x 3.0mm
9604361	Ultra II C18 3um	150 x 1.0mm
9604362	Ultra II C18 3um	150 x 2.1mm
9604362-700	Ultra II C18 3um	150 x 2.1mm With Trident
9604363	Ultra II C18 3um	150 x 3.2mm
9604363-700	Ultra II C18 3um	150 x 3.2mm With Trident
9604365	Ultra II C18 3um	150 x 4.6mm
9604365-700	Ultra II C18 3um	150 x 4.6mm With Trident
960436E	Ultra II C18 3um	150 x 3.0mm
96045	Ultra II C18 5um	Bulk Packing Material
96045	Ultra II C18 5um	Bulk Packing Material
96045	Ultra II C18 5um	Bulk Packing Material
960450210	Ultra II C18	10 x 4.0mm 3-pk
960450212	Ultra II C18	10 x 2.1mm 3-pk
960450220	Ultra II C18	20 x 4.0mm 2-pk
960450222	Ultra II C18	20 x 2.1mm 2-pk
9604510	Ultra II C18 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9604511	Ultra II C18 5um	100 x 1.0mm
9604512	Ultra II C18 5um	100 x 2.1mm

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034

ECHnology Pty Ltd

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

(Ctrl F)
> Cat.#s

> **Current System AUD Prices**

