

9604512-700	Ultra II C18 5um	100 x 2.1mm With Trident
9604513	Ultra II C18 5um	100 x 3.2mm
9604513-700	Ultra II C18 5um	100 x 3.2mm With Trident
9604515	Ultra II C18 5um	100 x 4.6mm
9604515-700	Ultra II C18 5um	100 x 4.6mm With Trident
9604517	Ultra II C18 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9604518	Ultra II C18 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9604519	Ultra II C18 5um	100 x 30mm Confirm customer has evaluated an analytical size column
960451E	Ultra II C18 5um	100 x 3.0mm
9604521	Ultra II C18 5um	200 x 1.0mm
9604522	Ultra II C18 5um	200 x 2.1mm
9604522-700	Ultra II C18 5um	200 x 2.1mm With Trident
9604523	Ultra II C18 5um	200 x 3.2mm
9604523-700	Ultra II C18 5um	200 x 3.2mm With Trident
9604525	Ultra II C18 5um	200 x 4.6mm
9604525-700	Ultra II C18 5um	200 x 4.6mm With Trident
960452E	Ultra II C18 5um	200 x 3.0mm
9604531	Ultra II C18 5um	30 x 1.0mm
9604532	Ultra II C18 5um	30 x 2.1mm
9604532-700	Ultra II C18 5um	30 x 2.1mm With Trident
9604533	Ultra II C18 5um	30 x 3.2mm
9604533-700	Ultra II C18 5um	30 x 3.2mm With Trident
9604535	Ultra II C18 5um	30 x 4.6mm
9604535-700	Ultra II C18 5um	30 x 4.6mm With Trident
960453E	Ultra II C18 5um	30 x 3.0mm
9604550	Ultra II C18 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9604551	Ultra II C18 5um	50 x 1.0mm
9604552	Ultra II C18 5um	50 x 2.1mm
9604552-700	Ultra II C18 5um	50 x 2.1mm With Trident
9604553	Ultra II C18 5um	50 x 3.2mm
9604553-700	Ultra II C18 5um	50 x 3.2mm With Trident
9604555	Ultra II C18 5um	50 x 4.6mm
9604555-700	Ultra II C18 5um	50 x 4.6mm With Trident
9604557	Ultra II C18 5um	50 x 10mm Confirm customer has evaluated an analytical size column
9604558	Ultra II C18 5um	50 x 21.2mm Confirm customer has evaluated an analytical size column
9604559	Ultra II C18 5um	50 x 30mm Confirm customer has evaluated an analytical size column
960455E	Ultra II C18 5um	50 x 3.0mm
9604560	Ultra II C18 5um	150 x 50mm Confirm customer has evaluated an analytical size column
9604561	Ultra II C18 5um	150 x 1.0mm
9604562	Ultra II C18 5um	150 x 2.1mm
9604562-700	Ultra II C18 5um	150 x 2.1mm With Trident
9604563	Ultra II C18 5um	150 x 3.2mm
9604563-700	Ultra II C18 5um	150 x 3.2mm With Trident
9604565	Ultra II C18 5um	150 x 4.6mm
9604565-700	Ultra II C18 5um	150 x 4.6mm With Trident
9604567	Ultra II C18 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9604568	Ultra II C18 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9604569	Ultra II C18 5um	150 x 30mm Confirm customer has evaluated an analytical size column
960456E	Ultra II C18 5um	150 x 3.0mm
9604570	Ultra II C18 5um	250 x 50mm Confirm customer has evaluated an analytical size column
9604571	Ultra II C18 5um	250 x 1.0mm
9604572	Ultra II C18 5um	250 x 2.1mm
9604572-700	Ultra II C18 5um	250 x 2.1mm With Trident
9604573	Ultra II C18 5um	250 x 3.2mm
9604573-700	Ultra II C18 5um	250 x 3.2mm With Trident
9604575	Ultra II C18 5um	250 x 4.6mm
9604575-700	Ultra II C18 5um	250 x 4.6mm With Trident
9604577	Ultra II C18 5um	250 x 10mm Confirm customer has evaluated an analytical size column
9604578	Ultra II C18 5um	250 x 21.2mm Confirm customer has evaluated an analytical size column
9604579	Ultra II C18 5um	250 x 30mm Confirm customer has evaluated an analytical size column
960457E	Ultra II C18 5um	250 x 3.0mm
9605311	Ultra II IBD 3um	100 x 1.0mm
9605312	Ultra II IBD 3um	100 x 2.1mm
9605312-700	Ultra II IBD 3um	100 x 2.1mm With Trident
9605313	Ultra II IBD 3um	100 x 3.2mm
9605313-700	Ultra II IBD 3um	100 x 3.2mm With Trident

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034

 <b>HROMalytic</b> +61(0)3 9762 2034 <b>ECHnology</b> Pty Ltd	<b>Australian Distributors</b> Importers & Manufacturers www.chromtech.net.au
Website NEW : <a href="http://www.chromalytic.net.au">www.chromalytic.net.au</a> E-mail : <a href="mailto:info@chromtech.net.au">info@chromtech.net.au</a> Tel: 03 9762 2034 . . . in AUSTRALIA	

[ Ctrl F ]  
> Cat.#s

> **Current System AUD Prices**

