

9605315	Ultra II IBD 3um	100 x 4.6mm
9605315-700	Ultra II IBD 3um	100 x 4.6mm With Trident
960531E	Ultra II IBD 3um	100 x 3.0mm
9605331	Ultra II IBD 3um	30 x 1.0mm
9605332	Ultra II IBD 3um	30 x 2.1mm
9605332-700	Ultra II IBD 3um	30 x 2.1mm With Trident
9605333	Ultra II IBD 3um	30 x 3.2mm
9605333-700	Ultra II IBD 3um	30 x 3.2mm With Trident
9605335	Ultra II IBD 3um	30 x 4.6mm
9605335-700	Ultra II IBD 3um	30 x 4.6mm With Trident
960533E	Ultra II IBD 3um	30 x 3.0mm
9605351	Ultra II IBD 3um	50 x 1.0mm
9605352	Ultra II IBD 3um	50 x 2.1mm
9605352-700	Ultra II IBD 3um	50 x 2.1mm With Trident
9605353	Ultra II IBD 3um	50 x 3.2mm
9605353-700	Ultra II IBD 3um	50 x 3.2mm With Trident
9605355	Ultra II IBD 3um	50 x 4.6mm
9605355-700	Ultra II IBD 3um	50 x 4.6mm With Trident
960535E	Ultra II IBD 3um	50 x 3.0mm
9605361	Ultra II IBD 3um	150 x 1.0mm
9605362	Ultra II IBD 3um	150 x 2.1mm
9605362-700	Ultra II IBD 3um	150 x 2.1mm With Trident
9605363	Ultra II IBD 3um	150 x 3.2mm
9605363-700	Ultra II IBD 3um	150 x 3.2mm With Trident
9605365	Ultra II IBD 3um	150 x 4.6mm
9605365-700	Ultra II IBD 3um	150 x 4.6mm With Trident
960536E	Ultra II IBD 3um	150 x 3.0mm
960550210	Ultra II IBD	10 x 4.0mm 3-pk
960550212	Ultra II IBD	10 x 2.1mm 3-pk
960550220	Ultra II IBD	20 x 4.0mm 2-pk
960550222	Ultra II IBD	20 x 2.1mm 2-pk
9605510	Ultra II IBD 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9605511	Ultra II IBD 5um	100 x 1.0mm
9605512	Ultra II IBD 5um	100 x 2.1mm
9605512-700	Ultra II IBD 5um	100 x 2.1mm With Trident
9605513	Ultra II IBD 5um	100 x 3.2mm
9605513-700	Ultra II IBD 5um	100 x 3.2mm With Trident
9605515	Ultra II IBD 5um	100 x 4.6mm
9605515-700	Ultra II IBD 5um	100 x 4.6mm With Trident
9605517	Ultra II IBD 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9605518	Ultra II IBD 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9605519	Ultra II IBD 5um	100 x 30mm Confirm customer has evaluated an analytical size column
960551E	Ultra II IBD 5um	100 x 3.0mm
9605521	Ultra II IBD 5um	200 x 1.0mm
9605522	Ultra II IBD 5um	200 x 2.1mm
9605522-700	Ultra II IBD 5um	200 x 2.1mm With Trident
9605523	Ultra II IBD 5um	200 x 3.2mm
9605523-700	Ultra II IBD 5um	200 x 3.2mm With Trident
9605525	Ultra II IBD 5um	200 x 4.6mm
9605525-700	Ultra II IBD 5um	200 x 4.6mm With Trident
960552E	Ultra II IBD 5um	200 x 3.0mm
9605531	Ultra II IBD 5um	30 x 1.0mm
9605532	Ultra II IBD 5um	30 x 2.1mm
9605532-700	Ultra II IBD 5um	30 x 2.1mm With Trident
9605533	Ultra II IBD 5um	30 x 3.2mm
9605533-700	Ultra II IBD 5um	30 x 3.2mm With Trident
9605535	Ultra II IBD 5um	30 x 4.6mm
9605535-700	Ultra II IBD 5um	30 x 4.6mm With Trident
960553E	Ultra II IBD 5um	30 x 3.0mm
9605550	Ultra II IBD 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9605551	Ultra II IBD 5um	50 x 1.0mm
9605552	Ultra II IBD 5um	50 x 2.1mm
9605552-700	Ultra II IBD 5um	50 x 2.1mm With Trident
9605553	Ultra II IBD 5um	50 x 3.2mm
9605553-700	Ultra II IBD 5um	50 x 3.2mm With Trident
9605555	Ultra II IBD 5um	50 x 4.6mm

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034

ECHnology Pty Ltd

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

(Ctrl F)
> Cat.#s

> **Current System AUD Prices**

