

960835E	Ultra II AQ C18 3um	50 x 3.0mm
9608361	Ultra II AQ C18 3um	150 x 1.0mm
9608362	Ultra II AQ C18 3um	150 x 2.1mm
9608362-700	Ultra II AQ C18 3um	150 x 2.1mm With Trident
9608363	Ultra II AQ C18 3um	150 x 3.2mm
9608363-700	Ultra II AQ C18 3um	150 x 3.2mm With Trident
9608365	Ultra II AQ C18 3um	150 x 4.6mm
9608365-700	Ultra II AQ C18 3um	150 x 4.6mm With Trident
960836E	Ultra II AQ C18 3um	150 x 3.0mm
960850210	Ultra II AQ C18	10 x 4.0mm 3-pk
960850212	Ultra II AQ C18	10 x 2.1mm 3-pk
960850220	Ultra II AQ C18	20 x 4.0mm 2-pk
960850222	Ultra II AQ C18	20 x 2.1mm 2-pk
9608510	Ultra II AQ C18 5um	100 x 50mm Confirm customer has evaluated an analytical size column
9608511	Ultra II AQ C18 5um	100 x 1.0mm
9608512	Ultra II AQ C18 5um	100 x 2.1mm
9608512-700	Ultra II AQ C18 5um	100 x 2.1mm With Trident
9608513	Ultra II AQ C18 5um	100 x 3.2mm
9608513-700	Ultra II AQ C18 5um	100 x 3.2mm With Trident
9608515	Ultra II AQ C18 5um	100 x 4.6mm
9608515-700	Ultra II AQ C18 5um	100 x 4.6mm With Trident
9608517	Ultra II AQ C18 5um	100 x 10mm Confirm customer has evaluated an analytical size column
9608518	Ultra II AQ C18 5um	100 x 21.2mm Confirm customer has evaluated an analytical size column
9608519	Ultra II AQ C18 5um	100 x 30mm Confirm customer has evaluated an analytical size column
960851E	Ultra II AQ C18 5um	100 x 3.0mm
9608521	Ultra II AQ C18 5um	200 x 1.0mm
9608522	Ultra II AQ C18 5um	200 x 2.1mm
9608522-700	Ultra II AQ C18 5um	200 x 2.1mm With Trident
9608523	Ultra II AQ C18 5um	200 x 3.2mm
9608523-700	Ultra II AQ C18 5um	200 x 3.2mm With Trident
9608525	Ultra II AQ C18 5um	200 x 4.6mm
9608525-700	Ultra II AQ C18 5um	200 x 4.6mm With Trident
960852E	Ultra II AQ C18 5um	200 x 3.0mm
9608531	Ultra II AQ C18 5um	30 x 1.0mm
9608532	Ultra II AQ C18 5um	30 x 2.1mm
9608532-700	Ultra II AQ C18 5um	30 x 2.1mm With Trident
9608533	Ultra II AQ C18 5um	30 x 3.2mm
9608533-700	Ultra II AQ C18 5um	30 x 3.2mm With Trident
9608535	Ultra II AQ C18 5um	30 x 4.6mm
9608535-700	Ultra II AQ C18 5um	30 x 4.6mm With Trident
960853E	Ultra II AQ C18 5um	30 x 3.0mm
9608550	Ultra II AQ C18 5um	50 x 50mm Confirm customer has evaluated an analytical size column
9608551	Ultra II AQ C18 5um	50 x 1.0mm
9608552	Ultra II AQ C18 5um	50 x 2.1mm
9608552-700	Ultra II AQ C18 5um	50 x 2.1mm With Trident
9608553	Ultra II AQ C18 5um	50 x 3.2mm
9608553-700	Ultra II AQ C18 5um	50 x 3.2mm With Trident
9608555	Ultra II AQ C18 5um	50 x 4.6mm
9608555-700	Ultra II AQ C18 5um	50 x 4.6mm With Trident
9608557	Ultra II AQ C18 5um	50 x 10mm Confirm customer has evaluated an analytical size column
9608558	Ultra II AQ C18 5um	50 x 21.2mm Confirm customer has evaluated an analytical size column
9608559	Ultra II AQ C18 5um	50 x 30mm Confirm customer has evaluated an analytical size column
960855E	Ultra II AQ C18 5um	50 x 3.0mm
9608560	Ultra II AQ C18 5um	150 x 50mm Confirm customer has evaluated an analytical size column
9608561	Ultra II AQ C18 5um	150 x 1.0mm
9608562	Ultra II AQ C18 5um	150 x 2.1mm
9608562-700	Ultra II AQ C18 5um	150 x 2.1mm With Trident
9608563	Ultra II AQ C18 5um	150 x 3.2mm
9608563-700	Ultra II AQ C18 5um	150 x 3.2mm With Trident
9608565	Ultra II AQ C18 5um	150 x 4.6mm
9608565-700	Ultra II AQ C18 5um	150 x 4.6mm With Trident
9608567	Ultra II AQ C18 5um	150 x 10mm Confirm customer has evaluated an analytical size column
9608568	Ultra II AQ C18 5um	150 x 21.2mm Confirm customer has evaluated an analytical size column
9608569	Ultra II AQ C18 5um	150 x 30mm Confirm customer has evaluated an analytical size column
960856E	Ultra II AQ C18 5um	150 x 3.0mm
9608570	Ultra II AQ C18 5um	250 x 50mm Confirm customer has evaluated an analytical size column

Chromalytic Technology Pty Ltd

AUSTRALIAN Distributors

03 9762 2034



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
 Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

(Ctrl F)
> Cat.#s

> Current System AUD Prices

