

Go To www.chromtech.net.au/filp/catalogs/"SUPPLIER" THEN add "/index.php" to URL
398 [Ctrl "F"] to SEARCH for Item > Description/Cat# etc

Typical Filtration PRODUCT LIST

prices as per June 2014
XR A11D/IISD-n 94

[Ctrl F]
> Cat.#s



> Current System AUD Prices

AUD

Norm-Ject Syringes Plastic Syringes All PP; Robust : LuerLok and Slip
in 1ml, 3ml, 5ml, 10ml, 20ml, 30,ml, 50ml

22766	Normject Plastic Syringes 1ml Tuberculin; Slip; 100-pk	48.00
22773	Normject Plastic Syringes 3ml; LuerLok; 100-pk	56.00
22775	Normject Plastic Syringes 5ml; LuerLok; 100-pk	88.00
22776	Normject Plastic Syringes 10ml; LuerLok; 100-pk	56.00
22777	Normject Plastic Syringes 30ml; LuerLok; 50-pk	78.00
22778	Normject Plastic Syringes 50ml; LuerLok; 30-pk	78.00

Syringe Filters

SF/MAT/Diam/PoreuM MAT = PTFE, PVDF, PolyProp, Nylon, PES, 25mmD ;100-pk (also 13mm)

SFPVFD025022N	0.22um (or 0.45um, non-sterile)	96.00
SFPVDF025022NG	Nonsterile PVDF Syringe Filters 0.22(um) 25(mm) GF prefilter	169.00
SFPVDF025045SM	Sterile G-MP Multi Layer PVDF 25mm 0.45um 1-/7/1/0.45um	195.00
6900-2502	Whatman GD/X : Sterile PVDF 0.22um multilayer. 25mm	420.00
6900-2504	Whatman GD/X : Sterile PVDF 0.45um multilayer, 25mm	420.00

Membrane Filters

MF/MAT/Diam/PoreuM MAT = PVDF, PolyProp, Nylon, PES; GlassFiber

MFPVDF047022N	0.22um or 0.45um; 47mmDiam; 200-pk	139.00
MFGF047100B	Glass Fiber Membrane Filter With Binder 1.00(um) 47(mm)	48.00
12-MV8010	Vacuum Pump : Hand Pump MityVAC	135.00
VPJ0332	Vacuum Pump /Motorised : Diaphragm 220VAC; 300mmHg 20L/Min	535.00

Filtration Systems : use separate MFs on top as needed

VFG030010	GlassFiltration System 47mm MF, Frit, Clamp, 300mlFunnel; 000ccBottle	295.00
	PES, PVDF (also MCE, CA, Nylon)	
VFPPVDF145250	Vac-uFill 47mmD, 250mlFunnel/250ml Bottle;10-pk	198.00
VFPPVDF222500	Vac-uFill 90mmD, 1.8umGF preFilter; 500mlFunnel/500ml Bottle, 10-pk	360.00

Cap Filtration

Vac-uCap, = Nalgene FAST Cap

CT-4622	VacuCAP;Pall; Filter 90mm/10pk Hydrophillic PES 0.2um; 10-pk	368.00
---------	--	--------

Liquid Pumps

NE-300AS	SyringePump Infusion Only NE-300AS	465.00
NE-1000AS	Syringe Pump Programmable NE-1000AS	825.00
Heater-Kit-1LG	Syringe Heater Kit ; PAD/controller & 120V 33W heating; inc transformer	545.00
NE-9000AS	Peristaltic Pump (Low Pressure Dispensing) NE-9000AS	725.00

Bottle Top Dispensers : also.5, 5.9. 10. 50, 50ml Max

BTDR-3	Microlit BottleTop Dispense 0.5to10ml Research Grade	215.00
--------	--	--------

Syringe Dispenser

CT-FDS-10	10ml SyringeDispenser Gun (Sterile)	45.00
C-4004	10ml Spring-loaded Syringe	18.00
DCV115-001	DoubleCheck Valve, Tubing Chimney Port RLL Exit CP2-5	3.95

Microfiltration : SingleStep Filtration; 200-pk Nylon; PES, PVDF, PTFE 0.2um, 0.45um

23964	0.45um PES	648.00
23960	0.45um PVDF	648.00

A:

**HROMalytic** +61(0)3 9762 2034
ECHnology Pty Ltd
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au