

Typical Vial & Accessories PRODUCT LIST


exGST prices as per
June 2014

[Ctrl F]
> Cat.#s

0.94
> Current System AUD Prices

Serum Pharma Vials CLEAR; 13mm or 20mm CrimpTop

CTSC-2V-13T-15D32H	2ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top	371	0.39
CTSC-4V-13T-16D38H	4ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top	425	0.39
CTSC-7V-20T-22D38H	7ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top	150	0.55
CTSC-10V-13T-22D45H	10ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top	198	0.66
CTSC-10V-20T-22D45H	10ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top *	324	0.56
CTSC-13V-20T-24D49H	13ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top	155	0.66
CTSC-13V-13T-24D49H	13ml Clear Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top	155	0.66
CTSC-20V-20T-26D66H	20ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top	145	0.66
CTSC-30V-20T-29D72H	30ml Clear Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top	113	0.77
CTSMCW-50V-20T-43D73H	50ml Clear Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top	113	0.85
CTSMCW-100V-20T-52D95H	50ml Clear Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top	81	0.95



Serum Pharma Vials AMBER; 13mm or 20mm CrimpTop

CTSA-2V-13T-15D32H	2ml Amber Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top	371	0.39
CTSA-3V-13T-15D35H	3ml Amber Neutral/Schott Borosilicate Type I Glass Vial 13mm Crimp Top	425	0.39
CTSA-7V-20T-22D38H	7ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top	210	0.55
CTSA-10V-20T-22D48H	10ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top	155	0.77
CTSA-13V-20T-24D50HTP	13ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top	155	0.77
CTSA-23V-20T-29D60HTP	23ml Amber Neutral/Schott Borosilicate Type I Glass Vial 20mm Crimp Top *	113	0.77
CTMAW-30V-20T-36D63H	30ml Amber Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top	92	0.89
CTMAW-50V-20T-43D73H	50ml Amber Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top	84	0.95
CTMAW-100V-20T-52D95H	100ml Amber Moulded Type II Glass Vial Sulphur Treated and Washed 20mm Crimp Top	66	0.95





Serum Stoppers; 13mm or 20mm Crimp Top ... also 32mm

CT-20BNRB-C	20mm Chlorobutyl Rubber Stopper	0.20
CT-20BNRB-FD2	20mm Grey Butyl Lyophilization Freeze Drying Rubber Stopper(2-Leg)	0.24
CT-13BNRB-C	13mm Chlorobutyl Rubber Stopper	0.20
CT-13BNRB-FD2	13mm Grey Butyl Lyophilization Freeze Drying Rubber Stopper(2-Leg)	0.24
CT-20BNRB	20mm Grey butyl rubber stopper (softer Light Gray type)	0.20
CT-20BNRB-R	20mm Butyl Rubber Stoppers - Red	0.20



CrimpTop Caps Plain Aluminium -X




CT-20CTO-Crimp-X	Aluminum Vial Seals 20mm Centre Tear-out : Color-X	0.24
CT-13CTO-Crimp-X	Aluminum Vial Seals 13mm Centre Tear-out : Color-X	0.24

Plastic FlipOff / WPFO -X "FLIP OFF" writing



Plastic FlipOff / KFO -X (plain Finish)



Hand Crimpers; Air, BenchPress, Hand Pneumatic Crimper with detachable Head ... also available



CT-Crimp-20S & -13S	HandCrimper 20mm (&13mm) Standard Plain AlumCaps	259.00
CT-Crimp-2OFO & -13FO	HandCrimper 20mm (&13mm) FlipOff Caps	259.00
CT-Decap-20 & -13	Hand Decapper 20mm (&13mm) both Plain and FO Cpas	259.00

Chromatography Vials 2ml 12x32ml : Autosampler; 8mm Universal : 9mm Agilent

ScrewCap LBSV012CSS LBSV012C LBSV032C LBSV012A LBSV032A LBSV002C LBSV022C LBSV002A LBSV022A LBSV102CSS CrimpTop LBSV042C LBSV062C LBSV042A LBSV062A LBSV112CRS	<p>Red PTFE/white silicone septa + Black screw cap with hole, /2ml 8-425 screw top vial,100-pk</p> <p>2ml Clear vial, 8-425 screw top,100-pk</p> <p>2ml Clear vial, 8-425 screw top, graduated with writing area,100-pk *</p> <p>2ml Amber vial, 8-425 screw top,100-pk</p> <p>2ml Amber vial, 8-425 screw top, graduated with writing area,100-pk</p> <p>2ml Clear vial, 9-425 screw top,100-pk</p> <p>2ml Clear vial, 9-425 screw top, graduated with writing area,100-pk *</p> <p>2ml Amber vial, 9-425 screw top,100-pk</p> <p>2ml Amber vial, 9-425 screw top, graduated with writing area,100-pk</p> <p>Red PTFE/white silicone septa + Blue screw cap with hole, /2ml 9-425 screw top vial,100-pk</p> <p>2ml Clear vial, crimp top,100-pk</p> <p>2ml Clear vial, crimp top, graduated with writing area,100-pk *</p> <p>2ml Amber vial, crimp top,100-pk</p> <p>2ml Amber vial, crimp top, graduated with writing area,100-pk</p> <p>White PTFE/red silicone septa+ Aluminium cap. for 2ml crimp top vial.100-pk</p>	 	<p>ASK for current prices</p>
--	---	------	---

Ask f

HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd
Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 ... in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au