



works with...BPAM (-001) Black Nylon over White Nylon

**Blood Pressure Monitor Bayonet Connector to Internal 5/16-24 UNF Thread and O-Ring**

PMSN (-XO) Panel Mount Nut 1/4-28 UNF with 7/16" Hex (For use with BPMB230)

**BLOOD PRESSURE:Dinamapr Monitor Adapters**

PM3N (-XO) Panel Mount Nut 5/16-24 UNF Zinc-Plated Class 2 Steel Nut (For use with BPMB240)

PM2N (-XO) Panel Mount Nut 3/8-24 UNF Zinc-Plated Class 2 Steel Nut (For Use with BPMB250)

**BLOOD PRESSURE:Bayonet Male Cuff Connectors**

BPM013 (-1) Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 1/8" (3.2 mm) ID Tubing

BPM240 (-1) Blood Pressure Monitor Bayonet Connector to 200 Series Barb, 5/32" (4.0 mm) ID Tubing

BPM025 (-1) Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 5/32" (4.0 mm) ID Tubing

BPM035 (-1) Blood Pressure Monitor Bayonet Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing

BPMMTL (-1 -5) Blood Pressure Monitor Bayonet Connector to Male Luer Integral Lock Ring

BPMB230 (-1) Blood Pressure Monitor Bayonet Panel Mount Connector to 200 Series Barb, 1/4-28 UNF Thread, 1/8" (3.2 mm) ID Tubing

**BLOOD PRESSURE:**

**Bayonet Panel Mount Connectors**

BPMB240 (-1) Blood Pressure Monitor Bayonet Panel Mount Connector to 200 Series Barb, 5/16-24 UNF Thread, 5/32" (4.0 mm) ID Tubing

BPMB250 (-1) Blood Pressure Monitor Bayonet Panel Mount Connector to 200 Series Barb, 3/8-24 UNF Thread, 3/16" (4.8 mm) ID Tubing

p47

**BLOOD PRESSURE**

**BLOOD PRESSURE: Dinamapr Cuff Connectors**

BPCL220 (-10 -81) Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Taper to 200 Series Barb, 3/32" (2.4 mm) ID Tubing

BPC220 (-10) Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 3/32" (2.4 mm) ID Tubing

BPCL230 (-10 -81) Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Internal Taper to 200 Series Barb, 1/8" (3.2 mm) ID Tubing

BPC230 (-10) Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 1/8" (3.2 mm) ID Tubing

BPCL025 (-10) Blood Pressure Adapter with 5/16-24 UNF Thread 6% Internal Taper to 500 Series Barb, 5/32" (4.0 mm) ID Tubing

BPC025 (-10) Blood Pressure Adapter with 5/16-24 UNF Thread to 500 Series Barb, 5/32" (4.0 mm) ID Tubing

BPCL240 (-10 -81) Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Internal Taper to 200 Series Barb, 5/32" (4.0 mm) ID Tubing

BPC240 (-10) Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 5/32" (4.0 mm) ID Tubing

BPCL250 (-10 -81) Blood Pressure Adapter with 5/16-24 UNF Thread and 6% Internal Taper to 200 Series Barb, 3/16" (4.8 mm) ID Tubing

BPC250 (-10) Blood Pressure Adapter with 5/16-24 UNF Thread to 200 Series Barb, 3/16" (4.8 mm) ID Tubing

CF229 (-1) Neonatal Blood Pressure Slip Fit (non luer) with 200 Series Barb, 3/32" (2.4 mm) ID Tubing

CF230 (-1) Neonatal Blood Pressure Slip Fit (non luer) with 200 Series Barb, 5/32" (4.0 mm) ID Tubing

p49

**TAPERED SEAL**

TSC2F007 (-6005) Female Taper Seal Non-Luer Connector to 500 Series Barb, 3/32" (2.4 mm) ID Tubing

TSC2F013 (-6005) Female Taper Seal Non-Luer Connector to 500 Series Barb, 1/8" (3.2 mm) ID Tubing

TSC2F035 (-6005) Female Taper Seal Non-Luer Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing

TSC2F210 (-6005) Female Taper Seal Non-Luer Connector to 200 Series Barb, 1/16" (1.6 mm) ID Tubing

MLB035 (-1 -8012 -9) Male Large Bore Non-Luer Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing

FLB035 (-1 -8012 -9) Female Large Bore Non-Luer Connector to 500 Series Barb, 3/16" (4.8 mm) ID Tubing

FLB055 (-1 -2 -8012 -9) Female Large Bore Non-Luer Connector to 500 Series Barb, 1/4" (6.4 mm) ID Tubing

FLB065 (-1 -8012 -9) Female Large Bore Non-Luer Connector to 500 Series Barb, 5/16" (8.0 mm) ID Tubing

FLBP (-1 -8012 -9) Female Large Bore Non-Luer Cap

MLBP (-1 -2 -8012 -9 -9002) Male Large Bore Non-Luer Plug

MLB055 (-1 -8012 -9)p Male Large Bore Non-Luer Connector to 500 Series Barb, 1/4" (6.4 mm) ID Tubing

MLB065 (-1 -8012 -9) Male Large Bore Non-Luer Connector to 500 Series Barb, 5/16" (8.0 mm) ID Tubing

MLBC7 (-9002) Male Large Bore Non-Luer Connector with Integral Lock Ring to Male Slip Luer (May be Used with Separate Rotating Lock Ring FSLLR)

MLBC8 (-9002) Male Large Bore Non-Luer Connector with Integral Lock Ring to Female Thread-Style Locking Luer

TSC2FB230 (-6005) Female Taper Seal Non-Luer with 1/4-28 UNF Panel Mount to 200 Series Barb, 1/8" (3.2 mm) ID Tubing

TSC2MP (-6005) 1,000 piece minimum

p51

**TUBE-TO-TUBE FITTINGS**

**TUBE-TO-TUBE: Straight Through**

N004 (-1 -2 -6005 -9 -J1A) Straight Through Connector with 500 Series Barbs, 1/16" (1.6 mm) ID Tubing

CC (-1 -2 -6005 -9 -J1A) Straight Through Connector with Classic Series Barbs, 1/8" (3.2 mm) ID Tubing

N025 (-1 -2 -6005 -9 -J1A) Straight Through Connector with 500 Series Barbs, 5/32" (4.0 mm) ID Tubing

N210 (-1 -2 -6005 -9 -J1A) Straight Through Connector with 200 Series Barbs, 1/16" (1.6 mm) ID Tubing

N410 (-1 -2 -6005 -J1A) Straight Through Connector with 400 Series Barbs, 1/16" (1.6 mm) ID Tubing

AA (-1 -2 -6005 -9 -J1A) Straight Through Connector with Classic Series Barbs, 1/16" (1.6 mm) ID Tubing

N007 (-1 -2 -6005 -9 -J1A) Straight Through Connector with 500 Series Barbs, 3/32" (2.4 mm) ID Tubing

N220 (-1 -2 -6005 -9 -9002 -J1A) Straight Through Connector with 200 Series Barbs, 3/32" (2.4 mm) ID Tubing

N420 (-1 -2 -6005 -J1A) Straight Through Connector with 400 Series Barbs, 3/32" (2.4 mm) ID Tubing

DD (-1 -2 -6005 -9 -J1A) Straight Through Connector with Classic Series Barbs, 3/32" (2.4 mm) ID Tubing