

BIOCHROM/HARVARD BIOSCIENCE/GE/PHARMACIA			
LIBRA S21, S22	L-990	P/XE	
ULTROSPEC 1000 E, I	LH-975	D2/W	INQUIRE
BIO-RAD			
BIOLOGIC DUAL FLOW			
760-1300	LH-25	D2	
760-1300	LHA-23	W	Lamp Ass'y
750-0216		HG	
1706, 1806A	L-295	D2	
1790, 1801, BIO-DIMENSION	LH-40	D2	
BIOFOCUS	LH-41	D2	
MICROSAMPLER 100	LH-41	D2	
REMEI	LH-115	D2	
SMARTSPEC 3000	L-60	P/XE	
BIOTAGE/DYAX/ARGONAUT			
HORIZON	L-735	HG	
ISOLERA SYSTEMS			
KILOPREP(100G)	LH-40	D2	
SP1, SP4, SPX	LH-60	D2	
FLASHMASTER II	LH-22	D2	
ISO-1SF, 1EF, 4SF, 4EF	L-735	HG	
ISO-1SV, 1SW, 1EV, 1EW	LH-60	D2	
ISO-4SV, 4SW, 4EV, 4EW	LH-60	D2	
BIOTEK INSTRUMENTS, INC			
EL-311	INQUIRE	W	
ELX-800 UV	L-1035	W	
ELX-808	INQUIRE	W	
FLX-800, 800T, 800TBI, 800TBID		W	
FLX-800TBX, 800 TC, 808/IU		W	
SYNERGY 2/AB, 2/FI/FP		W	
SYNERGY H4, HI/FI		W	
BISCHOFF			
LAMBDA 1000, 1010			INQUIRE
SP6, VISI-6			INQUIRE
BUCHI			
28062, 28063, 28067, 28068	LH-22	D2	
C-640			INQUIRE
BUCK			
BLC-20	LH-94	D2	
CAMAG			
TLC II SCANNER	L-285	D2	
CARLO ERBA/THERMO SCIENTIFIC			
	L-280	D2	
CARY/VARIAN			
1, 3, 4, 5, 100, 300, 400, 500	L-230	D2	
1, 3, 100, 300	L-930	W	
500	L-1015	W	
400	L-110	W	
50	L-1045	P/XE	
219	L-205	D2	
219	L-925	W	

CE 1200, 1220			(Ctrl F) > Cat.#s
CE 1020, 2020, 2040, 2041	L-115	D2	> Current System AUD Prices
CE 4000 thru 9000 SERIES	L-115	D2	
CE 2501, 2502, 7250	L-115	D2	
COLE -PALMER			
01953-00	LH-40	D2	
DIONEX/LC PACKING/GYNKOTEK			
AD20, PDA20	LH-5	D2	
AD25, PDA100, PDA2996	LH-3	D2	
PDA20		W	INQUIRE
PDA100		W	
CES I			INQUIRE
DSA-1	LH-40	D2	
FDM-2		P/XE	INQUIRE
RF1002/2000	L-800	XE	
UDM-170S, 340S			INQUIRE
UDM-37182			INQUIRE
UVD 160, 170S, 170U, 320, 340S, 340U	LH-91	D2	
VDM-2	L-235	D2	
DAD 3000, MWD 3000	LH-45	D2	
VWD 064377, MWD 064654	LH-45	D2	
VWD 3000	LH-45	D2	
D-STAR			
DDW-10, DLC-20, DVW-10	LH-94	D2	
ECOM			
LCD 2073A	LH-60	D2	
LCD 2071.2, LCD 2071.3	L-735	HG	
SAPPHIRE, OPAL, TOPAZ	LH-60	D2	
FLASH 06	LH-60	D2	
DYAX/SEE BIOTAGE			
EM SCIENCE/SEE HITACHI			
ERMA/ERC			
7511, 7512, 7515, 7522, 7530, 7535	L-610	W	
ESA			
520, 522	LH-40	D2	
FISHER /SEE THERMO SCIENTIFIC			
GBC			
AA	L-200	D2	
CINTRA 906, 908, 911	L-200	D2	
CINTRA 916		W	INQUIRE
CINTRA 1155, 1255	LH-910	P/XE	
CINTRA 1200	L-280	D2	
CINTRA 1205, 1206, 1210	LH-40	D2	
CINTRA 5100	LH-82	D2	