


GE HEALTHCARE/AMERSHAM/BIOSCIENCE/PHARMACIA/LKB			
2138, 2140, 2141			INQUIRE
2151	LH-20	D2	
4050, 4051, 4053, 4054	L-210	D2	
4050, 4060	L-835	W	
<b>AKTA SERIES</b>			
BASIC, EXPLORER, PILOT, PURIFIER	LH-990	P/XE	
AKTA MICROSYSTEM	LH-990	P/XE	
AKTA PRINME PLUS, PROCESS	LH-990	P/XE	
FPLC, AKTA PURIFIER 10, 100	L-1010 ***	HG	
ETAN LC 18-5050-50	LH-990	P/XE	
GENEQUANT II	LH-975	D2/W	INQUIRE
GENQUANT PRO			INQUIRE
RNA/DNA CALCULATOR			INQUIRE
NOVA SPEC II	L-835	W	
<b>ULTROSPEC SERIES</b>			
ULTROSPEC II, III, PLUS	L-210	D2	
ULTROSPEC II, III, PLUS	L-835	W	
ULTROSPEC 500 PRO	LH-125	D2/W	
ULTROSPEC 500 PRO	L-935	W	
ULTROSPEC 1000 e, 1000 i	LH-975	D2/W	INQUIRE
ULTROSPEC 2000, 3000, 4000	LH-125	D2/W	
ULTROSPEC 2000, 3000, 4000	L-935	W	
ULTROSPEC 3000 PRO	LH-125	D2/W	
ULTROSPEC 3000 PRO	L-935	W	
ULTROSPEC 1100 PRO	LH-975	D2/W	INQUIRE
ULTROSPEC 1100 PRO	L-935	W	
ULTROSPEC 2100 PRO, 3100 PRO	L-640	P/XE	
ULTROSPEC 3300, 4300	LH-125	D2/W	
ULTROSPEC 3300, 4300	L-935	W	
ULTROSPEC 5300 PRO, 6300 PRO	LH-125	D2/W	
UV-1	L-1030	HG	INQUIRE
UPC 900	L-1010 ***	HG	
UV 900	LH-990	P/XE	
***MUST RETURN SPENT LAMP AND HOUSING			
<b>GENTECH/SEE BEIJING RALEIGH</b>			
<b>GILFORD</b>			
200, 220, 222, 230, 240, 250, 252, 260	L-205	D2	
200, 220, 222, 230, 240, 250, 252, 260	L-875	W	
2000, 2400, 2600	L-205	D2	
2000, 2400, 2600	L-875	W	
FLUORO IV	L-801	XE	
RESPONSE I	L-210	D2	
RESPONSE I	L-900	W	
STASAR I, II, III		W	
SYVA S-3		W	
06/2014			

115, 116, 117			(Ctrl F) > Cat.#s	
118, 119, 151, 152, 155, 156			LH-151	D2
122			L-800	XE
<b>GRACE/ALLTECH</b>				
200, 201, 205			LH-40	D2
200, 201, 205			LH-760	W
305			LH-910	P/XE
REVELERIS UV1, UV2			LH-40	D2
<b>GROTON</b>				
FD100			L-890	P/XE
<b>GYNKOTEK/SEE DIONEX</b>				
<b>HACH</b>				
DR/4000U			LH-110	D2
DR/4000U, 4000V			L-730	W
DR/5000, DR/6000			LH-78	D2
DR/5000, DR/6000			L-780	W
48000, 48100			L-730	W
<b>HEWLETT-PACKARD/SEE AGILENT</b>				
<b>HITACHI</b>				
100-10, 20, 30, 40, 60, 80			L-280	D2
101, 102			L-430	D2
101, 102, 111, 124			L-805	W INQUIRE
124			L-210	D2
181, 191			L-805	W INQUIRE
150-20			LH-51	D2
150-20			L-805	W INQUIRE
200, 300, 330, 340			LH-51	D2
655A-22, 23			LH-51	D2
F-1000 THRU F-7000			L-800	XE
GENESPEC III			LH-58	D2
GENESPEC III				W INQUIRE
L-2400, 2420, 2450, 2455			LH-56	D2
L-2420, 2450, 2455			L-835	W
L-2480, 2485			L-800	XE
L-2500			LH-51	D2
L-3000, 3100, 3200, 4200, 4250			LH-51	D2
L-4000, 4500			LH-56	D2
L-4200, 4250			L-815	W
L-7000, 7300, 7400, 7420, 7450A, 7455			LH-56	D2
L-7420			L-815	W
L-7480, 7485, 7490			L-800	XE
SERIES 100			L-805	W INQUIRE
U-0080D			LH-56	D2
U-1100, 1500			LH-51	D2
U-1100			L-815	W