

MERCK/SEE HITACHI		
MICHROM BIORESOURCES/BRUKER		
UV-200	LH-40	D2
MD1-MD4	LH-25	D2
MICROTECH SCIENTIFIC, INC		
LINEAR MT200/200CE	LH-40	D2
LINEAR MT201, 205	LH-40	D2
MILTON ROY/SEE THERMO SCIENTIFIC/BAUSCH & LOMB		
MOLECULAR DEVICES		
SPECTRA MAX 340	L-65	P/XE
SPECTRA MAX 340 PC, 190, PLUS, 384	L-60	P/XE
SPECTRAMAX M2/M2E, M3, M4, M5/M5E	L-775	P/XE
SPECTRAMAX ANALYST GT	L-760	P/XE
SPECTRAMAX GEMINI , XPS, XS, EM	L-760	P/XE
SPECTRAMAX 250	L-760	P/XE
FLEX STATION, I, II, II 384	L-760	P/XE
NICOLET/SEE THERMO SCIENTIFIC/ATI		
OLIS		
RSM-1000	L-1000	XE
PERKIN-ELMER/EG&G		
HPLC UV DETECTORS		
785A	LH-98	D2
INTEGRAL 4000	L-210	D2
LC-55, 55B, 65, 75, 85B, 95, 135	L-210	D2
LC-200	LH-98	D2
LC-200		W INQUIRE
LC-235	L-210	D2
LC-480	LH-92	D2
SERIES 200/DAD, PHASE I & 2	LH-97	D2
SERIES 200/DAD	L-725	W
SERIES 200/FLUORO	L-800	XE
SERIES 200/UV-VIS/785A	LH-98	D2
FLUOROMETERS		
SERIES 200, 203, 204, 512, 650	L-800	XE
MPF-2A, 3, 4, 43, 43A, 44, 44A, 44B	L-800	XE
SPECTROPHOTOMETERS		
46, 55, 57	L-200	D2
220, 320, 330, 340	LH-51	D2
550, 551, 552, 555	L-240	D2
552A, 553, 559A	L-200	D2
570, 571, 572, 575, 576	L-200	D2
3840	L-200	D2
BIO	LH-92	D2
BIO	LH-845	W
EZ150, 201, 210	LH-51	D2
EZ150, 201, 210	L-840	W
06/2014		

PERKINELMER		<div> <div>> Cat.#s</div> <div> <div>Current System</div> <div>AUD Prices</div> </div> </div>	
LAMBDA 1, 3, 3B			
LAMBDA 1, 3B		LH-845	W
LAMBDA 2, 5, 7, 9, 10, 11, 12		LH-92	D2
LAMBDA 2, 5, 7, 9, 10, 11, 12		LH-845	W
LAMBDA 4, 6		LH-93	D2
LAMBDA 4, 6		LH-845	W
LAMBDA 14 THRU 19		LH-92	D2
LAMBDA 14 THRU 19		LH-845	W
LAMBDA 20, 25, 35, 40, 45		LH-92	D2
LAMBDA 20, 25, 35, 40, 45		LH-845	W
LAMBDA 650,750,800,850,900,950,1050		LH-92	D2
LAMBDA 650,750,800,850,900,950,1050		LH-845	W
ATOMIC ABSORPTION SPECTROMETERS			
100, 300 ANALYSIS		L-335	D2
200, 400 ANALYSIS			D2
272, 280, 360, 370, 372, 373, 380		L-310	D2
380, 460, 560		L-310	D2
410, 420, 430			D2
600, 700, 800 ANALYSIS		L-405	D2
1100, 1100B, 2100, 4100		L-405	D2
2280, 2380, 3030B		L-310	D2
3100, 3110, 3300 ANALYSIS		L-335	D2
4000 ANALYSIS		L-310	D2
5000, 5100		L-310	D2
PERSEPTIVE BIOSYSTEMS/SEE ABI			
PHARMACIA/SEE GE HEALTHCARE			
PHOTON TECHNOLOGY INT'L (PTI)			
QM-2/2005, QM-4/2005, QM-6/2005	L-970	XE	
QM-7/2005, QM-8/2005, QM-14/2005	L-970	XE	
POSTNOVA			
S3210, 3240	LH-13	D2	
PYE UNICAM/SEE THERMO SCIENTIFIC			
RAININ/SEE VARIAN			
ROCHE/APPLIED SCIENCE			
LIGHT CYCLE	L-1040	XE/R	
ROSEMONT ANALYTICAL			
890	L-330	D2	
SCHOEFFEL/ SEE ABI, PERKIN-ELMER			
SCIENTIFIC COMPUTING INTERNATIONAL			
FILMTEK			
4000, 4000-EM-DUV			
SCINCO			
S-3100, S-4100			
SEDEX			
55, 65, 75, 75C	L-935	W	
SEQUOIA TURNER			
240, 260, 260C		W	