

SHIMADZU			
HPLC UV DETECTORS			
LC2010, LC2010AHT, LC2010CHT	LH-84	D2	
SP-4	L-305	D2	
SPD-2A, 2AV, 6A, 6AV, 7A	L-305	D2	
SPD-6AV		W	INQUIRE
SPD-10A, 10AV, 10AVP, 10AVVP	LH-82	D2	
SPD-10AV, 10AVP, 10AVVP	L-830	W	
SPD-20A, 20AV	LH-82	D2	
SPD-20AV, SPD-M20AV	L-830	W	
SPD-M6A, M10A, M10AV		D2	INQUIRE
SPD-M6A, M10A,	L-935	W	
SPD-M10AV, M10AVP	L-245	W	
SPD-M10AVP	LH-88	D2	
SPD-M20A, M20AV	LH-87	D2	
SPD-M20A	L-245	W	
FLUOROMETERS			
RF-10A, 10AXL	L-800	XE	
RF-20A, 20AXS	L-795	XE	\$1350
RF- 535, 540, 551	L-800	XE	
RF-1501, 5000, 5301	L-800	XE	
SPECTROPHOTOMETERS			
BIOSPEC 1601	L-265	D2	
BIOSPEC 1601	L-110	W	
CS-920, 930,,	L-265	D2	
MULTISPEC 1501	L-265	D2	
MULTISPEC 1501	L-935	W	
UV-120, 120-01, 120-02, 140-02	L-265	D2	
UV-120	L-110	W	
UV-120-01, 120-02,		W	INQUIRE
UV-150, 150-02, 160A, 160U, 163	L-265	D2	
UV-160A, 160U	L-110	W	
UV-240, 250, 260, 265,	L-265	D2	
UV-240, 260, 265	L-110	W	
UV-1200, 1201, 1201PC, 1203	L-265	D2	
UV-1200, 1201, 1203	L-245	W	
UV-1240, 1501, 1600, 1601, 1601PC	L-265	D2	
UV-1240, 1240V, 1700, 1800	L-245	W	
UV-1501, 1600, 1601, 1601PC	L-110	W	
UV-1603, 1650, 1700, 1800, 2100	L-265	D2	
UV-1603, 1650	L-110	W	
2101, 2401PC, 2450, 2600, 2700	L-265	D2	
UV-2101PC, 2401, 2450, 2600, 2700	L-110	W	
UV-2501PC, 2550, 3100, 3100PC	L-265	D2	
UV-2501PC, 2550,	L-110	W	
UV-3101PC, 3600, 3700	L-265	D2	
UV-3101PC, 3600	L-110	W	
ATOMIC ABSORPTION SPECTROMETERS			
AA-625, 630, 640, 660, 670, 680U	L-265	D2	
AA-6200	L-265	D2	

HPLC Detector & Spectrometer UV Lamps 2014

06/2014

RI 73, 74, 101			
SIEVERS/IONICS INSTR. BUS.GROUP/GE ANALYT .INST.			
400, 800	L-700	HG	
900, 5310C, PPT	L-605	HG	
SKALAR			
6750	LH-22		
SPECTRA-PHYSICS/SEE THERMO SCIENTIFIC			
SPECTRONICS/SEE THERMO SCIENTIFIC			
SPEX/SEE JY (JOBIN YVON)			
TECAN			
SPECTRO FLUOR PLUS	L-1025	P/XE	\$1095
GENIOS	L-1025	P/XE	\$1095
INFINITE F200, M200 NANQUANT	L-1025	P/XE	\$1095
INFINITE F200 PRO, M200 PRO	L-1025	P/XE	\$1095
THERMO SCIENTIFIC			
AMINCO			
AB-2 (SERIES 2)			INQUIRE
SLM 2000	L-280	D2	
SLM 3000	L-270	D2	
SLM 4800, 8000, 8100	L-855	XE	
ATI			
ATI/240, 250	LHA-34	Lamp Ass'y	
	LH-34	Cartridge Only	
BAUSCH & LOMB			
AQUAMATE, AQUAMATE UV	L-910	W	
AQUAMATE UV	LH-54	D2	
BIOMATE 3		P/XE	INQUIRE
BIOMATE 5	L-910	W	
DE 50	L-225	D2	
EPSILON	L-915	W	
GENESYS 2, 5		D2	INQUIRE
GENESYS 2, 5	L-915	W	
GENESYS 6, 60, 10 BIO, 10 UV		P/XE	INQUIRE
GENESYS 8	LH-54	D2	
GENESYS 8	L-910	W	
GENESYS 10, 20	L-915	W	
EVOLUTION 60		P/XE	INQUIRE
EVOLUTION 60 S			INQUIRE
SPECTRONIC 20			INQUIRE
SPECTRONIC 21, -DV, MV, UVD, UVM	L-945	W	
SPECTRONIC 21D, 21 DUV	L-225	D2	
SPECTRONIC 200, 601, 700, 701, 710	L-225	D2	
SPECTRONIC 70, 100, 700, 701, 710	L-920	W	
SPECTRONIC 301, 401	L-915	W	
SPECTRONIC 501, 601 ROUND TOP	L-615	W	
SPECTRONIC 501, 601 POINTED TOP	L-915	W	
SPECTRONIC 1001, 1201, 2000	L-225	D2	
SPECTRONIC 1001, 1201	L-810	W	
SPECTRONIC 2000	L-920	W	
ENVIRONMENTAL INSTRUMENTS/SO2 ANALYZER			
43A, 43CT, 43CH	L-40	P/XE	
FISHER			
20-0344	LH-40	D2	